

ABX Minocal



(Exp.) 2023-07-05

(YYYY - MM - DD)

LOT		CX 483		CAL		(Exp.) 2023-07-05 (YYYY - MM - DD)						
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA		
		ES60	CRP 200			60 8P 60 18P						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	9.0			9.0					± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.47	4.49			4.47					± 0.06	
HB HGB	g/dl	13.2	13.0			13.2					± 0.2	
	g/l	132	130			132					± 2	
	mmol/l	8.20	8.09			8.20					± 0.12	
HT HCT	%	37.1	37.3			37.8					± 1.0	
	l/l	0.371	0.373			0.378					± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	229	247			243					± 10	
VMP MPV	µm ³ ; fl	8.6	8.3			8.8					± 0.5	

ABX Minocal



LOT		CX 483		CAL		(Exp.) 2023-07-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
		MICROS	MICROS								
		Scil Vet ABC	Scil Vet ABC Plus ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.8								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.49	4.49								± 0.06
HB HGB	g/dl	13.1	13.1								± 0.2
	g/l	131	131								± 2
	mmol/l	8.14	8.14								± 0.12
HT HCT	%	37.3	37.5								± 1.0
	l/l	0.373	0.375								± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	229	243								± 10
VMP MPV	µm ³ ; fl	N/A	8.6								± 0.5

Ref: TEMP-0387 Rev.49 FRONT / RECTO 1300113555