

ABX Minocal

LOT

CX 483

CAL



(Exp.) 2023-07-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	9.0			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.47	4.49			4.47					± 0.06
HB HGB	g/dl	13.2	13.0			13.2					± 0.2
	g/l	132	130			132					± 2
	mmol/l	8.20	8.09			8.20					± 0.12
HT HCT	%	37.1	37.3			37.8					± 1.0
	l/l	0.371	0.373			0.378					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	229	247			243					± 10
VMP MPV	µm ³ ; fl	8.6	8.3			8.8					± 0.5

ABX Minocal

LOT

CX 483

CAL



(Exp.) 2023-07-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.8	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.49	4.49							± 0.06
HB HGB	g/dl	13.1	13.1							± 0.2
	g/l	131	131							± 2
	mmol/l	8.14	8.14							± 0.12
HT HCT	%	37.3	37.5							± 1.0
	l/l	0.373	0.375							± 0.010
PLA PLT	10 ⁹ /mm ³ ;10 ⁹ /l	229	243							± 10
VMP MPV	µm ³ ; fl	N/A	8.6							± 0.5

Ref: TEMP-0387 Rev.49 FRONT / RECTO 1300113555