

ABX Minocal



(Exp.) 2023-07-05

(YYYY - MM - DD)

LOT		CX 483		CAL		(Exp.) 2023-07-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.0			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.49			4.52					± 0.06
HB HGB	g/dl	13.3	13.3			13.3					± 0.2
	g/l	133	133			133					± 2
	mmol/l	8.26	8.26			8.26					± 0.12
HT HCT	%	37.7	37.3			37.7					± 1.0
	l/l	0.377	0.373			0.377					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	251	247			245					± 10
VMP MPV	µm ³ ; fl	8.2	8.3			7.4					± 0.5

Ref. TEMP-0387 Rev.49 FRONT / RECTO 1300090308