

ABX Minocal



LOT

CX 483

CAL



(Exp.) 2023-07-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3				
		H550	H500 OT H500 CT			H500 OT	H500 CT H550			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.90	8.90			8.90	8.90			± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.52			4.50	4.50			± 0.06
HB HGB	g/dl	12.8	12.8			12.8	12.8			± 0.2
	g/l	128	128			128	128			± 2
	mmol/l	7.95	7.95			7.95	7.95			± 0.12
HT HCT	%	37.1	37.1			35.5	34.8			± 1.0
	l/l	0.371	0.371			0.355	0.348			± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	258	258			283	288			± 10
VMP MPV	µm ³ ; fl	10.3	10.3			N/A	N/A			± 0.5

Ref: TEMP-0387 Rev.49 BACK / VERSO 1300113555