

ABX Minocal



(Exp.) 2023-08-05

(YYYY - MM - DD)

LOT		CX 484		CAL		(Exp.) 2023-08-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.9			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.52			4.53					± 0.06
HB HGB	g/dl	13.7	13.6			13.7					± 0.2
	g/l	137	136			137					± 2
	mmol/l	8.51	8.45			8.51					± 0.12
HT HCT	%	39.0	38.2			39.0					± 1.0
	l/l	0.390	0.382			0.390					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	248	247			246					± 10
VMP MPV	µm ³ ; fl	8.5	8.4			7.9					± 0.5

 1300090308
 FRONT / RECTO
 Ref. TEMP-0387 Rev.49