

ABX Minocal



| LOT                      |                                                       | CX 484      |         | CAL    |        | <div><div></div><div>(Exp.) 2023-08-05</div><div>(YYYY - MM - DD)</div></div> |        |        |        |        |                         |
|--------------------------|-------------------------------------------------------|-------------|---------|--------|--------|-------------------------------------------------------------------------------|--------|--------|--------|--------|-------------------------|
|                          |                                                       | Rev 1       |         |        |        |                                                                               |        |        |        |        |                         |
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | Minilysebio |         |        |        |                                                                               |        |        |        |        | TOLERANCES<br>TOLERANCE |
|                          |                                                       | MICROS      | MICROS  | MICROS | MICROS | MICROS                                                                        | PENTRA | PENTRA | PENTRA | PENTRA |                         |
|                          |                                                       | ES60        | CRP 200 |        |        | 60 8P<br>60 18P                                                               |        |        |        |        |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 9.0         | 8.9     |        |        | 8.9                                                                           |        |        |        |        | ± 0.2                   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.53        | 4.52    |        |        | 4.51                                                                          |        |        |        |        | ± 0.06                  |
| HB HGB                   | g/dl                                                  | 13.8        | 13.3    |        |        | 13.6                                                                          |        |        |        |        | ± 0.2                   |
|                          | g/l                                                   | 138         | 133     |        |        | 136                                                                           |        |        |        |        | ± 2                     |
|                          | mmol/l                                                | 8.57        | 8.28    |        |        | 8.45                                                                          |        |        |        |        | ± 0.12                  |
| HT HCT                   | %                                                     | 38.7        | 38.2    |        |        | 39.0                                                                          |        |        |        |        | ± 1.0                   |
|                          | l/l                                                   | 0.387       | 0.382   |        |        | 0.390                                                                         |        |        |        |        | ± 0.010                 |
| PLA PLT                  | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 236         | 247     |        |        | 244                                                                           |        |        |        |        | ± 10                    |
| VMP MPV                  | µm <sup>3</sup> ; fl                                  | 9.0         | 8.4     |        |        | 9.3                                                                           |        |        |        |        | ± 0.5                   |

# ABX Minocal



| LOT                      |                                                       | CX 484       |                             | CAL |  | <div><div></div><div>(Exp.) 2023-08-05</div><div>(YYYY - MM - DD)</div></div> |  |  |  |  |                         |
|--------------------------|-------------------------------------------------------|--------------|-----------------------------|-----|--|------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--|--|--|-------------------------|
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | LyseVet      |                             |     |  |                                                                                                                                                                  |  |  |  |  | TOLERANCES<br>TOLERANCE |
|                          |                                                       | MICROS       | MICROS                      |     |  |                                                                                                                                                                  |  |  |  |  |                         |
|                          |                                                       | Scil Vet ABC | Scil Vet ABC Plus<br>ESV 60 |     |  |                                                                                                                                                                  |  |  |  |  |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 8.8          | 8.8                         |     |  |                                                                                                                                                                  |  |  |  |  | ± 0.2                   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.55         | 4.52                        |     |  |                                                                                                                                                                  |  |  |  |  | ± 0.06                  |
| HB HGB                   | g/dl                                                  | 13.5         | 13.3                        |     |  |                                                                                                                                                                  |  |  |  |  | ± 0.2                   |
|                          | g/l                                                   | 135          | 133                         |     |  |                                                                                                                                                                  |  |  |  |  | ± 2                     |
|                          | mmol/l                                                | 8.38         | 8.26                        |     |  |                                                                                                                                                                  |  |  |  |  | ± 0.12                  |
| HT HCT                   | %                                                     | 38.9         | 38.4                        |     |  |                                                                                                                                                                  |  |  |  |  | ± 1.0                   |
|                          | l/l                                                   | 0.389        | 0.384                       |     |  |                                                                                                                                                                  |  |  |  |  | ± 0.010                 |
| PLA PLT                  | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 234          | 257                         |     |  |                                                                                                                                                                  |  |  |  |  | ± 10                    |
| VMP MPV                  | µm <sup>3</sup> ; fl                                  | N/A          | 9.0                         |     |  |                                                                                                                                                                  |  |  |  |  | ± 0.5                   |

Ref: TEMP-0387 Rev.49 FRONT / RECTO 1300113555