

ABX Minocal

LOT

CX 484

CAL



(Exp.) 2023-08-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.9			8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.52			4.51					± 0.06
HB HGB	g/dl	13.8	13.3			13.6					± 0.2
	g/l	138	133			136					± 2
	mmol/l	8.57	8.28			8.45					± 0.12
HT HCT	%	38.7	38.2			39.0					± 1.0
	l/l	0.387	0.382			0.390					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	236	247			244					± 10
VMP MPV	µm ³ ; fl	9.0	8.4			9.3					± 0.5

ABX Minocal

LOT

CX 484

CAL



(Exp.) 2023-08-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.8	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.55	4.52							± 0.06
HB HGB	g/dl	13.5	13.3							± 0.2
	g/l	135	133							± 2
	mmol/l	8.38	8.26							± 0.12
HT HCT	%	38.9	38.4							± 1.0
	l/l	0.389	0.384							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	234	257							± 10
VMP MPV	µm ³ ; fl	N/A	9.0							± 0.5

Ref: TEMP-0387 Rev.49 FRONT / RECTO 1300113555