

ABX Erytrol



LOT

NX 442

Rev 1

CONTROL



(Exp.) 2023-09-05

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL		1	TOLERANCES TOLERANCE	CONTROL		2	TOLERANCES TOLERANCE	CONTROL		3	TOLERANCES TOLERANCE
		PENTRA	PENTRA	PENTRA		PENTRA	PENTRA	PENTRA					
		DX120	DX NEXUS	DX120	DX NEXUS	DX120	DX NEXUS						
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.5	7.5		± 0.8	18.3	18.3		± 1.5	2.6	2.6		± 0.4
GBC CWBC	10 ³ /mm ³ ; 10 ⁹ /l	7.3	7.3		± 0.8	16.8	16.8		± 1.4	1.6	1.6		± 0.4
ERB	%	2.8	2.8		± 1.7	9.5	9.5		± 3.3	67.5	67.5		± 23.6
ERB	#	0.204	0.204		± 0.100	1.596	1.596		± 0.700	1.080	1.080		± 0.378

1219930076-C

Ref: TEMP-0822 Rev.36

GBC / CWBC

* les valeurs sont données à titre d'indication et doivent être utilisées comme RuO (Utilisation en Recherche Seulement) Indicative expected values, must be used as RuO (Research used only)