

ABX Minocal


 (Exp.) **2023-08-05**
 (YYYY - MM - DD)

LOT		CX 484				CAL		(Exp.) 2023-08-05 (YYYY - MM - DD)			
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60	DX120 - DF120	60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	9.0	9.1	8.9	9.2	9.0	9.2	9.1	8.9	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48	4.47	4.48	4.53	4.56	4.49	4.56	4.54	4.53	± 0.06
HB HGB	g/dl	13.4	13.4	13.5	13.6	13.6	13.5	13.6	13.6	13.6	± 0.2
	g/l	134	134	135	136	136	135	136	136	136	± 2
	mmol/l	8.32	8.32	8.38	8.45	8.45	8.38	8.45	8.45	8.45	± 0.12
HT HCT	%	38.1	38.0	37.0	37.6	37.4	36.8	37.4	37.5	37.6	± 1.0
	l/l	0.381	0.380	0.370	0.376	0.374	0.368	0.374	0.375	0.376	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	249	249	234	243	253	245	253	249	243	± 10
VMP MPV	µm ³ ; fl	8.8	8.8	N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.49 BACK / VERSO 1300090308