

ABX Minocal



LOT

CX 484

CAL



(Exp.) 2023-08-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3				
		H550	H500 OT H500 CT			H500 OT	H500 CT H550			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80			8.80	8.80			± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.57			4.54	4.54			± 0.06
HB HGB	g/dl	13.2	13.2			13.3	13.3			± 0.2
	g/l	132	132			133	133			± 2
	mmol/l	8.20	8.20			8.26	8.26			± 0.12
HT HCT	%	38.6	38.6			36.9	36.4			± 1.0
	l/l	0.386	0.386			0.369	0.364			± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	253	253			272	277			± 10
VMP MPV	μm ³ ; fl	10.4	10.4			N/A	N/A			± 0.5

Ref: TEMP-0387 Rev.49 BACK / VERSO 1300113555