



Yumizen G PT Liq 4 for Yumizen G Line



Reference No.:		1300036373		
Lot No.:		930411 (rev. 1)		
Expiry date:		2024-10-31		
Coagulometers		ISI	MASTER CURVE	
Optical	Yumizen G Line		1,18	%
		100		12,8
		50		17,9
		33		23,1
		25		28,0
	MNPT	12,8		
Important notice				
<p>The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.</p>				