



LAtrol™ Normal Control Plasma ***LAtrol®***

HORIBA Reference No.:	1300079435
Biomedica Reference No.:	816N
Lot No.:	220803
Expiry date:	2024/08/03

for Yumizen G800

Parameter	Unit	Range
DRVVT	sec	39.2 – 52.2
DRVVC	sec	34.1 – 43.1
Normalized Ratio	-	<1.16

for Yumizen G1500/G1550

Parameter	Unit	Range
DRVVT	sec	-
DRVVC	sec	-
Normalized Ratio	-	<1.16

Important notice

Each laboratory should establish its own DRVVT and DRVVC normal reference range (mean value ± 2 S.D.).

The Normalized Ratio is calculated as follows:

$$\text{Normalized Ratio} = \frac{\text{DRVVT (sample)} / \text{DRVVT (mean clot time test)}}{\text{DRVVC (sample)} / \text{DRVVC (mean clot time confirm)}}$$