

# ABX Minocal



(Exp.) 2023-09-05

(YYYY - MM - DD)

| LOT                      |                                                       | CX 485          |         | CAL    |        | <div><div></div><div>(Exp.) 2023-09-05</div><div>(YYYY - MM - DD )</div></div> |        |        |        |        |                         |
|--------------------------|-------------------------------------------------------|-----------------|---------|--------|--------|--------------------------------------------------------------------------------|--------|--------|--------|--------|-------------------------|
|                          |                                                       | Rev 1           |         |        |        |                                                                                |        |        |        |        |                         |
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | ABX Lyse        |         |        |        |                                                                                |        |        |        |        | TOLERANCES<br>TOLERANCE |
|                          |                                                       | MICROS          | MICROS  | MICROS | MICROS | MICROS                                                                         | PENTRA | PENTRA | PENTRA | PENTRA |                         |
|                          |                                                       | ES60<br>Care ST | CRP 200 |        |        | 60 8P<br>60 18P<br>ADVIA 60                                                    |        |        |        |        |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 9.1             | 9.2     |        |        | 9.3                                                                            |        |        |        |        | ± 0.2                   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.55            | 4.61    |        |        | 4.57                                                                           |        |        |        |        | ± 0.06                  |
| HB HGB                   | g/dl                                                  | 13.5            | 13.6    |        |        | 13.4                                                                           |        |        |        |        | ± 0.2                   |
|                          | g/l                                                   | 135             | 136     |        |        | 134                                                                            |        |        |        |        | ± 2                     |
|                          | mmol/l                                                | 8.38            | 8.45    |        |        | 8.32                                                                           |        |        |        |        | ± 0.12                  |
| HT HCT                   | %                                                     | 38.4            | 38.5    |        |        | 38.2                                                                           |        |        |        |        | ± 1.0                   |
|                          | l/l                                                   | 0.384           | 0.385   |        |        | 0.382                                                                          |        |        |        |        | ± 0.010                 |
| PLA PLT                  | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 240             | 236     |        |        | 239                                                                            |        |        |        |        | ± 10                    |
| VMP MPV                  | µm <sup>3</sup> ; fl                                  | 8.3             | 8.2     |        |        | 7.4                                                                            |        |        |        |        | ± 0.5                   |

Ref: TEMP-0387 Rev.49 FRONT / RECTO 1300090308