

ABX Minocal



LOT		CX 485		CAL		 (Exp.) 2023-09-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.2			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.61			4.58					± 0.06
HB HGB	g/dl	13.5	13.3			13.6					± 0.2
	g/l	135	133			136					± 2
	mmol/l	8.38	8.28			8.45					± 0.12
HT HCT	%	38.5	38.5			38.5					± 1.0
	l/l	0.385	0.385			0.385					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	233	236			238					± 10
VMP MPV	µm ³ ; fl	8.5	8.2			8.6					± 0.5

ABX Minocal



LOT		CX 485		CAL		(Exp.) 2023-09-05 (YYYY - MM - DD)					
		Rev 1	LyseVet								
		MICROS	MICROS								
PARAMETRES PARAMETERS		UNITES UNITS	Scil Vet ABC	Scil Vet ABC Plus ESV 60							TOLERANCES TOLERANCE
GB WBC		10 ³ /mm ³ ;10 ⁹ /l	8.9	9.0							± 0.2
GR RBC		10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.55							± 0.06
HB HGB		g/dl	13.3	13.3							± 0.2
		g/l	133	133							± 2
		mmol/l	8.26	8.26							± 0.12
HT HCT		%	38.2	38.2							± 1.0
		l/l	0.382	0.382							± 0.010
PLA PLT		10 ³ /mm ³ ;10 ⁹ /l	216	232							± 10
VMP MPV		µm ³ ; fl	N/A	8.4							± 0.5

Ref: TEMP-0387 Rev.49 FRONT / RECTO 1300113555