

ABX Minocal



LOT

CX 485

CAL



(Exp.) 2023-09-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3				
		H550	H500 OT H500 CT			H500 OT	H500 CT H550			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.10	9.10			9.10	9.10			± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.58			4.53	4.53			± 0.06
HB HGB	g/dl	13.0	13.0			13.0	13.0			± 0.2
	g/l	130	130			130	130			± 2
	mmol/l	8.07	8.07			8.07	8.07			± 0.12
HT HCT	%	38.0	38.0			36.6	36.1			± 1.0
	l/l	0.380	0.380			0.366	0.361			± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	246	246			275	270			± 10
VMP MPV	μm ³ ; fl	10.0	10.0			N/A	N/A			± 0.5

Ref: TEMP-0387 Rev.49 BACK / VERSO 1300113555