

# Yumizen G PT 5 for Yumizen G Line



Reference No.:		1300036338		
Lot No.:		930405 (rev. 1)		
Expiry date:		2025-04-30		
Coagulometers		ISI	MASTER CURVE	
Optical	Yumizen G Line		%	sec
		100	12,8	
		50	17,7	
		33	22,7	
		25	27,5	
		1,18	MNPT	12,8
<b>Important notice</b>				
<ul style="list-style-type: none"> <li>- The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.</li> <li>- The efficiency of the reconstitution of the reagent may be improved by performing it at 37°C.</li> </ul>				