



# ***Yumizen G FIB 2 for Yumizen G1500/G1550 and Yumizen G800***



<b>Reference No.:</b>		<b>1300036383</b>	
<b>Lot No.:</b>		<b>930315 (rev. 1)</b>	
<b>Expiry date:</b>		<b>2025-03-31</b>	
<b>Coagulometers</b>		<b>MASTER CURVE</b>	
<b>Optical</b>	Yumizen G1500/G1550/G800	<b>g/L</b>	<b>sec</b>
		<b>4</b>	<b>9,0</b>
		<b>3</b>	<b>10,9</b>
		<b>2</b>	<b>14,2</b>
		<b>1</b>	<b>22,2</b>
	<b>Dilution</b>	<b>10x</b>	

# ***Yumizen G FIB 2 for Yumizen G400 DDi/G200***



<b>Reference No.:</b>		<b>1300036383</b>	
<b>Lot No.:</b>		<b>930315 (rev. 1)</b>	
<b>Expiry date:</b>		<b>2025-03-31</b>	
<b>Coagulometers</b>		<b>MASTER CURVE</b>	
<b>Optical</b>	Yumizen G400 DDi/G200	<b>g/L</b>	<b>sec</b>
		<b>4</b>	<b>9,7</b>
		<b>3</b>	<b>11,7</b>
		<b>2</b>	<b>15,4</b>
		<b>1</b>	<b>24,4</b>
	<b>Dilution</b>	<b>10x</b>	