

ABX Minocal

LOT

CX 486

CAL



(Exp.) 2023-10-05

Rev **2**

(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.0			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.51			4.53					± 0.06
HB HGB	g/dl	13.3	13.0			13.4					± 0.2
	g/l	133	130			134					± 2
	mmol/l	8.26	8.09			8.32					± 0.12
HT HCT	%	37.4	37.4			37.8					± 1.0
	l/l	0.374	0.374			0.378					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	227	240			245					± 10
VMP MPV	µm ³ ; fl	8.2	8.1			8.5					± 0.5

ABX Minocal

LOT

CX 486

Rev 2

CAL



(Exp.) **2023-10-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.49							± 0.06
HB HGB	g/dl	13.2	13.0							± 0.2
	g/l	132	130							± 2
	mmol/l	8.20	8.07							± 0.12
HT HCT	%	37.5	37.3							± 1.0
	l/l	0.375	0.373							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	225	239							± 10
VMP MPV	µm ³ ; fl	N/A	8.5							± 0.5

Ref: TEMP-0387 Rev.50 FRONT / RECTO 1300113555