

ABX Minocal

LOT
CX 486

 Rev **2**
CAL

(Exp.) 2023-10-05

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.0			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.51			4.53					± 0.06
HB HGB	g/dl	13.3	13.0			13.4					± 0.2
	g/l	133	130			134					± 2
	mmol/l	8.26	8.09			8.32					± 0.12
HT HCT	%	37.4	37.4			37.8					± 1.0
	l/l	0.374	0.374			0.378					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	227	240			245					± 10
VMP MPV	µm ³ ; fl	8.2	8.1			8.5					± 0.5

ABX Minocal



LOT		CX 486		CAL								 (Exp.) 2023-10-05 (YYYY - MM - DD)		
		Rev 2		LyseVet										
PARAMETRES PARAMETERS	UNITES UNITS											TOLERANCES TOLERANCE		
		MICROS		MICROS										
		Scil Vet ABC	Scil Vet ABC Plus ESV 60											
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.7									± 0.2		
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.49									± 0.06		
HB HGB	g/dl	13.2	13.0									± 0.2		
	g/l	132	130									± 2		
	mmol/l	8.20	8.07									± 0.12		
HT HCT	%	37.5	37.3									± 1.0		
	l/l	0.375	0.373									± 0.010		
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	225	239									± 10		
VMP MPV	µm ³ ; fl	N/A	8.5									± 0.5		

Ref: TEMP-0387 Rev.50 FRONT / RECTO 1300113555