

# ABX Minocal

**LOT**

**CX 486**

**CAL**



**(Exp.) 2023-10-05**

(YYYY-MM-DD)

Rev **2**

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.0</b>	<b>9.0</b>			<b>9.1</b>					± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.51</b>	<b>4.51</b>			<b>4.53</b>					± 0.06
HB HGB	g/dl	<b>13.3</b>	<b>13.0</b>			<b>13.4</b>					± 0.2
	g/l	<b>133</b>	<b>130</b>			<b>134</b>					± 2
	mmol/l	<b>8.26</b>	<b>8.09</b>			<b>8.32</b>					± 0.12
HT HCT	%	<b>37.4</b>	<b>37.4</b>			<b>37.8</b>					± 1.0
	l/l	<b>0.374</b>	<b>0.374</b>			<b>0.378</b>					± 0.010
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>227</b>	<b>240</b>			<b>245</b>					± 10
VMP MPV	µm <sup>3</sup> ; fl	<b>8.2</b>	<b>8.1</b>			<b>8.5</b>					± 0.5

# ABX Minocal

**LOT**

**CX 486**

Rev 2

**CAL**



(Exp.) **2023-10-05**  
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.0</b>	<b>8.7</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.52</b>	<b>4.49</b>							± 0.06
HB HGB	g/dl	<b>13.2</b>	<b>13.0</b>							± 0.2
	g/l	<b>132</b>	<b>130</b>							± 2
	mmol/l	<b>8.20</b>	<b>8.07</b>							± 0.12
HT HCT	%	<b>37.5</b>	<b>37.3</b>							± 1.0
	l/l	<b>0.375</b>	<b>0.373</b>							± 0.010
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>225</b>	<b>239</b>							± 10
VMP MPV	µm <sup>3</sup> ; fl	N/A	<b>8.5</b>							± 0.5

Ref: TEMP-0387 Rev.50 FRONT / RECTO 1300113555