

ABX Minocal


 (Exp.) 2023-10-05
 (YYYY - MM - DD)

LOT		CX 486		CAL		(Exp.) 2023-10-05 (YYYY - MM - DD)						
		Rev 2										
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA		
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60	DX120 - DF120	60 60 C+ ES60	XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.1	9.2	9.1	9.3	9.1	9.3	9.1	9.1	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.51	4.45	4.51	4.51	4.48	4.51	4.50	4.51	± 0.06	
HB HGB	g/dl	13.1	13.1	13.2	13.3	13.2	13.1	13.2	13.2	13.3	± 0.2	
	g/l	131	131	132	133	132	131	132	132	133	± 2	
	mmol/l	8.14	8.14	8.20	8.26	8.20	8.14	8.20	8.20	8.26	± 0.12	
HT HCT	%	37.4	37.0	36.0	36.5	36.1	35.8	36.1	36.2	36.5	± 1.0	
	l/l	0.374	0.370	0.360	0.365	0.361	0.358	0.361	0.362	0.365	± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	237	242	235	243	245	241	245	241	243	± 10	
VMP MPV	µm ³ ; fl	8.3	8.3	N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5	

Ref: TEMP-0387 Rev.50 BACK / VERSO 1300090308