

LOT		MX 443 L		Rev 1		CONTROL		L		 (Exp.)		2023-12-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	2.1	2.0			2.1						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.33	2.32			2.29						± 0.17	
HB HGB	g/dl	6.1	6.1			6.1						± 0.5	
	g/l	61	61			61						± 5	
	mmol/l	3.79	3.79			3.79						± 0.31	
HT HCT	%	16.8	16.5			16.3						± 2.0	
	l/l	0.168	0.165			0.163						± 0.020	
VGM MCV	µm ³ ; fl	72	71			71						± 4	
TGMH MCH	pg	26.2	26.3			26.6						± 2.0	
	fmol	1.63	1.63			1.65						± 0.12	
CCMH MCHC	g/dl	36.4	37.0			37.5						± 3.0	
	g/l	364	370			375						± 30	
	mmol/l	22.58	23.00			23.30						± 1.86	
IDR RDW	%	15.0	14.0			14.0						± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	77	75			76						± 22	
VPM MPV	µm ³ ; fl	8.2	8.0			7.5						± 2.0	
LYMPHO	#	1.33	1.28			1.35						± 0.40	
	%	64.5	65.0			65.5						± 8.0	
MONO	#	0.20	0.20			0.20						± 0.20	
	%	8.5	9.0			8.5						± 6.0	
GRANULO	#	0.57	0.52			0.55						± 0.40	
	%	27.0	26.0			26.0						± 7.0	

LOT		MX 443 N		CONTROL		N		 (Exp.)		2023-12-05 (YYYY - MM - DD)											
		Rev 1																			
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse								TOLERANCES TOLERANCE									
				MICROS		MICROS		MICROS		MICROS				ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP 200								60 8P - 60 18P ADVIA 60							
GB	WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.6	7.3					7.5									± 0.8			
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.52	4.48					4.49									± 0.26			
HB	HGB	g/dl	13.5	13.4					13.5									± 0.7			
		g/l	135	134					135								± 7				
		mmol/l	8.38	8.32					8.38								± 0.43				
HT	HCT	%	37.1	36.7				36.8									± 2.7				
		l/l	0.371	0.367				0.368									± 0.027				
VGM	MCV	µm ³ ; fl	82	82				82									± 4				
TGMH	MCH	pg	29.9	29.9				30.1									± 2.0				
		fmol	1.85	1.86				1.87									± 0.12				
CCMH	MCHC	g/dl	36.4	36.5				36.7									± 3.0				
		g/l	364	365				367									± 30				
		mmol/l	22.62	22.65				22.77									± 1.86				
IDR	RDW	%	14.5	13.5				13.5									± 3.0				
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l	251	251				246									± 40				
VPM	MPV	µm ³ ; fl	8.0	7.5				7.5									± 2.0				
LYMPHO	#	2.51	2.41				2.48										± 0.40				
	%	33.0	33.0				33.0										± 5.0				
MONO	#	0.57	0.55				0.53										± 0.30				
	%	7.5	7.5				7.0										± 4.0				
GRANULO	#	4.52	4.34				4.50										± 0.60				
	%	59.5	59.5				60.0										± 6.0				

Ref: TEMP-0830 Rev.42 FRONT / RECTO 1300104144

LOT		MX 443 H		CONTROL		H		 (Exp.)		2023-12-05 (YYYY - MM - DD)											
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse								TOLERANCES TOLERANCE									
				MICROS		MICROS		MICROS		MICROS				ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP 200				60 8P - 60 18P ADVIA 60											
GB	WBC	10 ³ /mm ³ ; 10 ⁹ /l		20.1		20.2				20.2										± 1.7	
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l		5.70		5.66				5.69										± 0.32	
HB	HGB	g/dl		18.0		18.0				18.0										± 0.9	
		g/l		180		180				180								± 9			
		mmol/l		11.18		11.18				11.18								± 0.56			
HT	HCT	%		50.2		49.8				50.6										± 3.4	
		l/l		0.502		0.498				0.506										± 0.034	
VGM	MCV	µm ³ ; fl		88		88				89										± 4	
TGMH	MCH	pg		31.6		31.8				31.6										± 2.0	
		fmol		1.96		1.97				1.96										± 0.12	
CCMH	MCHC	g/dl		35.9		36.1				35.5										± 3.0	
		g/l		359		361				355										± 30	
		mmol/l		22.28		22.44				22.07										± 1.86	
IDR	RDW	%		14.0		13.0				12.5										± 3.0	
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l		492		481				488										± 65	
VPM	MPV	µm ³ ; fl		7.7		7.5				7.0										± 2.0	
LYMPHO	#			3.52		3.43				3.33										± 1.00	
	%			17.5		17.0				16.5										± 5.0	
MONO	#			1.01		1.21				1.01										± 0.80	
	%			5.0		6.0				5.0										± 4.0	
GRANULO	#			15.58		15.55				15.86										± 1.40	
	%			77.5		77.0				78.5										± 7.0	