

ABX Minocal

LOT		CX 487		CAL		(Exp.) 2023-11-05 (YYYY - MM - DD)						
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA		
		ES60	CRP 200			60 8P 60 18P						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	8.9			9.3					± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.56			4.54					± 0.06	
HB HGB	g/dl	13.4	13.1			13.5					± 0.2	
	g/l	134	131			135					± 2	
	mmol/l	8.32	8.15			8.38					± 0.12	
HT HCT	%	37.6	38.1			38.1					± 1.0	
	l/l	0.376	0.381			0.381					± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	220	230			237					± 10	
VMP MPV	µm ³ ; fl	8.8	8.4			8.9					± 0.5	

ABX Minocal



LOT		CX 487		CAL		(Exp.) 2023-11-05 (YYYY - MM - DD)							
		Rev 1		LyseVet									
PARAMETRES PARAMETERS	UNITES UNITS										TOLERANCES TOLERANCE		
		MICROS		MICROS									
		Scil Vet ABC	Scil Vet ABC Plus ESV 60										
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.9								± 0.2		
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.53								± 0.06		
HB HGB	g/dl	13.2	13.1								± 0.2		
	g/l	132	131								± 2		
	mmol/l	8.20	8.14								± 0.12		
HT HCT	%	37.9	37.6								± 1.0		
	l/l	0.379	0.376								± 0.010		
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	216	237								± 10		
VMP MPV	µm ³ ; fl	N/A	8.8								± 0.5		

Ref: TEMP-0387 Rev.50 FRONT / RECTO 1300113555