

ABX Minocal

LOT

CX 487

CAL



(Exp.) 2023-11-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	8.9			9.3					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.56			4.54					± 0.06
HB HGB	g/dl	13.4	13.1			13.5					± 0.2
	g/l	134	131			135					± 2
	mmol/l	8.32	8.15			8.38					± 0.12
HT HCT	%	37.6	38.1			38.1					± 1.0
	l/l	0.376	0.381			0.381					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	220	230			237					± 10
VMP MPV	µm ³ ; fl	8.8	8.4			8.9					± 0.5

ABX Minocal

LOT

CX 487

Rev 1

CAL



(Exp.) **2023-11-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.53							± 0.06
HB HGB	g/dl	13.2	13.1							± 0.2
	g/l	132	131							± 2
	mmol/l	8.20	8.14							± 0.12
HT HCT	%	37.9	37.6							± 1.0
	l/l	0.379	0.376							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	216	237							± 10
VMP MPV	µm ³ ; fl	N/A	8.8							± 0.5

Ref: TEMP-0387 Rev.50 FRONT / RECTO 1300113555