

ABX Minocal


 (Exp.) 2023-11-05
 (YYYY - MM - DD)

LOT		CX 487		CAL		(Exp.) 2023-11-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60	DX120 - DF120	60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.2	9.3	9.3	9.4	9.2	9.4	9.2	9.3	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.53	4.46	4.52	4.56	4.51	4.56	4.55	4.52	± 0.06
HB HGB	g/dl	13.2	13.2	13.2	13.3	13.3	13.2	13.3	13.3	13.3	± 0.2
	g/l	132	132	132	133	133	132	133	133	133	± 2
	mmol/l	8.20	8.20	8.20	8.26	8.26	8.20	8.26	8.26	8.26	± 0.12
HT HCT	%	37.2	37.1	36.1	36.6	36.3	35.9	36.3	36.6	36.6	± 1.0
	l/l	0.372	0.371	0.361	0.366	0.363	0.359	0.363	0.366	0.366	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	233	241	225	226	237	229	237	231	226	± 10
VMP MPV	µm ³ ; fl	8.8	8.9	N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.50 BACK / VERSO 1300090308