

ABX Minocal

LOT

CX 488

CAL



(Exp.) 2023-12-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	8.9			9.3					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.35			4.42					± 0.06
HB HGB	g/dl	13.1	12.8			13.1					± 0.2
	g/l	131	128			131					± 2
	mmol/l	8.14	7.97			8.14					± 0.12
HT HCT	%	36.7	36.8			37.3					± 1.0
	l/l	0.367	0.368			0.373					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	227	240			236					± 10
VMP MPV	µm ³ ; fl	8.9	8.7			9.0					± 0.5

ABX Minocal

LOT

CX 488

Rev 1

CAL



(Exp.) **2023-12-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.40	4.39							± 0.06
HB HGB	g/dl	13.0	12.9							± 0.2
	g/l	130	129							± 2
	mmol/l	8.07	8.01							± 0.12
HT HCT	%	37.4	37.1							± 1.0
	l/l	0.374	0.371							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	219	231							± 10
VMP MPV	µm ³ ; fl	N/A	9.0							± 0.5

Ref: TEMP-0387 Rev.50 FRONT / RECTO 1300113555