

ABX Erytrol



LOT

NX 444

Rev 1

CONTROL



(Exp.) 2024-01-05

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL		1	TOLERANCES TOLERANCE	CONTROL		2	TOLERANCES TOLERANCE	CONTROL		3	TOLERANCES TOLERANCE
		PENTRA	PENTRA	DX120		DX NEXUS	DX120	DX NEXUS		DX120	DX NEXUS		
		DX120	DX NEXUS										
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.4	7.4		± 0.8	17.7	17.7		± 1.5	2.7	2.7		± 0.4
GBC CWBC	10 ³ /mm ³ ; 10 ⁹ /l	7.1	7.1		± 0.8	16.0	16.0		± 1.4	1.5	1.5		± 0.4
ERB	%	4.3	4.3		± 1.7	9.9	9.9		± 3.3	79.3	79.3		± 50.0
ERB	#	0.305	0.305		± 0.100	1.584	1.584		± 0.700	1.190	1.190		± 0.500

1219930076-C

Ref: TEMP-0822 Rev.36

GBC / CWBC

* les valeurs sont données à titre d'indication et doivent être utilisées comme RuO (Utilisation en Recherche Seulement) Indicative expected values, must be used as RuO (Research used only)