

ABX Minocal



(Exp.) 2023-12-05

(YYYY - MM - DD)

LOT		CX 488		CAL								
		Rev 1				ABX Lyse						
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS		MICROS		PENTRA		PENTRA		TOLERANCES TOLERANCE
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	8.9			9.1						± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.35			4.38						± 0.06
HB HGB	g/dl	13.2	13.1			13.1						± 0.2
	g/l	132	131			131						± 2
	mmol/l	8.20	8.14			8.14						± 0.12
HT HCT	%	37.1	36.8			37.0						± 1.0
	l/l	0.371	0.368			0.370						± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	233	240			233						± 10
VMP MPV	µm ³ ; fl	8.7	8.7			8.0						± 0.5

Ref: TEMP-0387 Rev.50 FRONT / RECTO 1300090308