

ABX Minocal


 (Exp.) 2023-12-05
 (YYYY - MM - DD)

LOT		CX 488		CAL		(Exp.) 2023-12-05 (YYYY - MM - DD)					
Rev 1		ABX Lysebio									
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60	DX120 - DF120	60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.2	9.5	9.2	9.4	9.0	9.4	9.2	9.2	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.37	4.39	4.38	4.39	4.39	4.34	4.39	4.40	4.39	± 0.06
HB HGB	g/dl	12.9	13.0	13.0	13.1	13.0	13.0	13.0	13.0	13.1	± 0.2
	g/l	129	130	130	131	130	130	130	130	131	± 2
	mmol/l	8.01	8.07	8.07	8.14	8.07	8.07	8.07	8.07	8.14	± 0.12
HT HCT	%	36.3	36.7	35.5	35.8	35.3	35.2	35.3	35.4	35.8	± 1.0
	l/l	0.363	0.367	0.355	0.358	0.353	0.352	0.353	0.354	0.358	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	235	241	230	235	240	233	240	235	235	± 10
VMP MPV	µm ³ ; fl	8.9	8.8	N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.50 BACK / VERSO 1300090308