

ABX Minotrol Retic

LOT		RX 444		RETIC ADJUSTMENT TOOL				(Exp.) 2024-01-05 (YYYY - MM - DD)			
Rev 1											
PARAMETERS	UNITS										
		LEVEL			1	TOLERANCE	LEVEL			3	TOLERANCE
		YUMIZEN			YUMIZEN						
		H2500 Since V1.7			H2500 Since V1.7						
RET	%	1.70			± 0.40					N/A	
MFI	-				N/A	0.479			± 0.030		