


# ABX Minotrol Retic

**LOT** RX 444  
Rev 1

**RETIC ADJUSTMENT TOOL**

 (Exp.) **2024-01-05**  
(YYYY - MM - DD)

PARAMETERS	UNITS	LEVEL 1				LEVEL 3			
		YUMIZEN				YUMIZEN			
		H2500 Since V1.7			TOLERANCE	H2500 Since V1.7			TOLERANCE
RET	%	<b>1.70</b>			± 0.40				N/A
MFI	-				N/A	<b>0.479</b>			± 0.030