

ABX Minocal



LOT

CX 488

CAL



(Exp.) 2023-12-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE	
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3					
		H550	H500 OT H500 CT			H500 OT	H500 CT H550				
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.95	8.95			8.98	8.98				± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.41	4.41			4.39	4.39				± 0.06
HB HGB	g/dl	12.7	12.7			12.7	12.7				± 0.2
	g/l	127	127			127	127				± 2
	mmol/l	7.89	7.89			7.89	7.89				± 0.12
HT HCT	%	37.0	37.0			35.8	35.2				± 1.0
	l/l	0.370	0.370			0.358	0.352				± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	239	239			250	261				± 10
VMP MPV	μm ³ ; fl	10.2	10.2			N/A	N/A				± 0.5

Ref: TEMP-0387 Rev.50 BACK / VERSO 1300113555