

ABX Minocal



LOT		CX 489		CAL		(Exp.) 2024-01-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.6			8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48	4.53			4.43					± 0.06
HB HGB	g/dl	13.3	13.1			13.3					± 0.2
	g/l	133	131			133					± 2
	mmol/l	8.26	8.15			8.26					± 0.12
HT HCT	%	37.6	38.3			37.9					± 1.0
	l/l	0.376	0.383			0.379					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	220	235			238					± 10
VMP MPV	µm ³ ; fl	8.6	8.4			9.0					± 0.5

ABX Minocal



LOT		CX 489		CAL		(Exp.)2024-01-05 (YYYY - MM - DD)								
		Rev 1		LyseVet										
PARAMETRES PARAMETERS	UNITES UNITS											TOLERANCES TOLERANCE		
		MICROS		MICROS										
		Scil Vet ABC	Scil Vet ABC Plus ESV 60											
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.6									± 0.2		
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.44	4.48									± 0.06		
HB HGB	g/dl	13.1	13.1									± 0.2		
	g/l	131	131									± 2		
	mmol/l	8.14	8.14									± 0.12		
HT HCT	%	38.0	37.9									± 1.0		
	l/l	0.380	0.379									± 0.010		
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	219	226									± 10		
VMP MPV	µm ³ ; fl	N/A	8.9									± 0.5		

Ref: TEMP-0387 Rev.50 FRONT / RECTO 1300113555