

ABX Minocal

LOT

CX 489

CAL



(Exp.) 2024-01-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.6			8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48	4.53			4.43					± 0.06
HB HGB	g/dl	13.3	13.1			13.3					± 0.2
	g/l	133	131			133					± 2
	mmol/l	8.26	8.15			8.26					± 0.12
HT HCT	%	37.6	38.3			37.9					± 1.0
	l/l	0.376	0.383			0.379					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	220	235			238					± 10
VMP MPV	µm ³ ; fl	8.6	8.4			9.0					± 0.5

ABX Minocal

LOT

CX 489

CAL



(Exp.) 2024-01-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.6							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.44	4.48							± 0.06
HB HGB	g/dl	13.1	13.1							± 0.2
	g/l	131	131							± 2
	mmol/l	8.14	8.14							± 0.12
HT HCT	%	38.0	37.9							± 1.0
	l/l	0.380	0.379							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	219	226							± 10
VMP MPV	µm ³ ; fl	N/A	8.9							± 0.5

Ref: TEMP-0387 Rev.50 FRONT / RECTO 1300113555