

ABX Minocal



LOT		CX 489		CAL		(Exp.) 2024-01-05 (YYYY - MM - DD)					
Rev 1		ABX Lyse									
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	8.6			8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.45	4.53			4.45					± 0.06
HB HGB	g/dl	13.3	13.4			13.4					± 0.2
	g/l	133	134			134					± 2
	mmol/l	8.26	8.32			8.32					± 0.12
HT HCT	%	38.0	38.3			37.8					± 1.0
	l/l	0.380	0.383			0.378					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	231	235			230					± 10
VMP MPV	µm ³ ; fl	8.5	8.4			7.8					± 0.5

Ref: TEMP-0387 Rev.50 FRONT / RECTO 1300090308

 1300090308
 Ref. TEMP-0387 Rev.50 FRONT / RECTO