

ABX Minocal



LOT		CX 489		CAL				(Exp.) 2024-01-05 (YYYY - MM - DD)					
		Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE		
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA			
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60	DX120 - DF120	60 60 C+ ES60	XL80			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.0	9.1	8.9	9.2	8.9	9.2	9.0	8.9	± 0.2		
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.46	4.43	4.40	4.45	4.47	4.41	4.47	4.49	4.45	± 0.06		
HB HGB	g/dl	13.2	13.2	13.2	13.3	13.3	13.2	13.3	13.3	13.3	± 0.2		
	g/l	132	132	132	133	133	132	133	133	133	± 2		
	mmol/l	8.20	8.20	8.20	8.26	8.26	8.20	8.26	8.26	8.26	± 0.12		
HT HCT	%	37.7	37.4	36.3	36.7	36.4	35.9	36.4	36.6	36.7	± 1.0		
	l/l	0.377	0.374	0.363	0.367	0.364	0.359	0.364	0.366	0.367	± 0.010		
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	237	238	229	228	239	229	239	232	228	± 10		
VMP MPV	µm ³ ; fl	9.0	8.7	N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5		

1300090308
BACK / VERSO
Ref. TEMP-0387 Rev.50