

ABX Minocal



LOT

CX 489

CAL



(Exp.) 2024-01-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE	
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3					
		H550	H500 OT H500 CT			H500 OT	H500 CT H550				
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80			8.80	8.80				± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.50			4.45	4.45				± 0.06
HB HGB	g/dl	13.0	13.0			13.1	13.1				± 0.2
	g/l	130	130			131	131				± 2
	mmol/l	8.07	8.07			8.14	8.14				± 0.12
HT HCT	%	37.8	37.8			36.5	36.4				± 1.0
	l/l	0.378	0.378			0.365	0.364				± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	247	247			252	259				± 10
VMP MPV	µm ³ ; fl	10.7	10.7			N/A	N/A				± 0.5

Ref: TEMP-0387 Rev.50 BACK / VERSO 1300113555