

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 1 of 21

Homocysteine R1

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Product Name: Homocysteine R1

Product code: H7575-R1

1.2 Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses: For in vitro quantitative determination of total homocysteine in human serum and plasma.

Uses advised against: Not determined or not applicable.

Reasons why uses advised against: Not determined or not applicable.

1.3 Details of the manufacturer/supplier of the safety data sheet

Manufacturer:

United States

HORIBA Instruments Incorporated

5449 Research Drive

Canton, MI 48188

734-487-8300

horiba.com

1.4 Emergency telephone number:

United States

HORIBA Instruments Incorporated

1-800-445-9853 (24 hours per day)

France

Organisme de conseil/centre antipoison national

+33 1 45 42 59 59 (24 hours per day)

Portugal

Órgão consultor nacional/Centro Antivenenos

+351 800 250 250 (24 hours per day)

Spain

Centro de información toxicológica/organismo asesor nacional

+34 91 562 04 20 (24 hours per day)

Czech Republic

Národní poradní orgán/toxikologické středisko

+420 224 919 293 (24 hours per day)

Greece

Εθνικό συμβουλευτικό όργανο/Κέντρο Δηλητηριάσεων

+30 210 779 3777 (24 hours per day)

Romania

Organism consultativ național/Centru pentru otrăviri

+40 21 3183606 (24 hours per day)

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 2 of 21

Homocysteine R1

Poland

Krajowa instytucja doradcza/Ośrodek zatruć
+48 22 619 66 54 (24 hours per day)

SECTION 2: Hazard(s) identification

2.1 Classification of the substance or mixture:

Classification according to Regulation (EC) No. 1272/2008 (CLP):

Acute toxicity (oral), category 4

Hazard-determining components of labeling:

2-amino-2-(hydroxymethyl)propane-1,3-diol

Sodium azide

Ethanol

Additional Information: None

2.2 Label elements

Labelling according to Regulation (EC) No 1272/2008 (CLP)

Hazard pictograms:



Signal Word: Warning

Hazard statements:

H302 Harmful if swallowed

Precautionary statements:

P264 Wash skin thoroughly after handling.

P270 Do not eat, drink or smoke when using this product

P301+P312 IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.

P330 Rinse mouth

P501 Dispose of contents into sewer system after diluting with large volumes of water, if in accordance with local regulations.

2.3 Other hazards: None known

SECTION 3: Composition/information on ingredients

3.1 Substance: Not applicable.

3.2 Mixture:

Identification	EU REACH Registration No.	Name	Classification according to Regulation (EC) No. 1272/2008 (CLP)	Weight %
CAS number: 77-86-1 EC number: 201-064-4	-	2-amino-2-(hydroxymethyl)propane-1,3-diol	Skin Irrit. 2; H315 STOT SE 3 (RI); H335 Eye Irrit. 2; H319	1-10
CAS number: 64-17-5 EC number: 200-578-6	-	Ethanol	Flam. Liq. 2; H225 Eye Irrit. 2; H319	1-10

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 3 of 21

Homocysteine R1

CAS number: 26628-22-8 EC number: 247-852-1	-	Sodium azide	Acute Tox. 2 (Oral); H300 Aquatic Acute 1; H400 Acute Tox. 1 (Dermal); H310 Acute Tox. 2 (Inh); H330 Aquatic Chronic 1; H410 STOT RE 2; H373 M-Factor: 1 EUH032	<1
--	---	--------------	--	----

Additional information: None

Full Text of H and EUH statements: See section 16

SECTION 4: First aid measures

4.1 Description of first aid measures

General notes:

Show this Safety Data Sheet to the doctor in attendance. Take precautions to ensure your own safety before attempting rescue. Wear appropriate safety eyewear, gloves, protective clothing and respiratory protection to prevent exposure. See Section 8 of this SDS for personal protective equipment recommendations. Do not use the mouth to mouth method if victim has ingested or inhaled the product. Give artificial respiration with the aid of a pocket mask equipped with a one-way valve or other proper device.

Following inhalation:

If inhaled, remove person to fresh air and place in a position comfortable for breathing. Keep person at rest. If breathing is difficult, administer oxygen. If breathing has stopped, provide artificial respiration. If symptoms develop or persist, seek medical advice/attention.

Following skin contact:

Remove contaminated clothing and shoes. Rinse skin with copious amounts of water [shower] for several minutes. Launder contaminated clothing before reuse. If symptoms develop or persist, seek medical advice/attention.

Following eye contact:

Rinse eyes with plenty of water for several minutes. Remove contact lenses, if present and easy to do so. Protect unexposed eye. If symptoms develop or persist, seek medical advice/attention.

Following ingestion:

If swallowed, DO NOT induce vomiting unless told to do so by a physician or poison control center. Rinse mouth with water. Never give anything by mouth to an unconscious person. If spontaneous vomiting occurs, place on the left side with head down to prevent aspiration of liquid into the lungs. If symptoms develop or persist, seek medical advice/attention.

Self-Protection of the first aider:

Not determined or not available.

4.2 Most important symptoms and effects, both acute and delayed

Acute symptoms and effects:

Acute oral exposure may lead to dizziness, drowsiness, headache, breathing difficulties, nausea, vomiting, abdominal pain, and lowering of consciousness. Adverse effects are dependent on exposure (dose, concentration, contact time).

Delayed symptoms and effects:

Symptoms of exposure may be delayed.

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 4 of 21

Homocysteine R1

4.3 Indication of any immediate medical attention and special treatment needed

Specific treatment:

Not determined or not available.

Notes for the doctor:

Treat symptomatically.

SECTION 5: Firefighting measures

5.1 Extinguishing media

Suitable extinguishing media:

Water mist/fog, carbon dioxide, dry chemical or alcohol resistant foam.

Unsuitable extinguishing media:

Do not use water jet.

5.2 Special hazards arising from the substance or mixture:

Thermal decomposition may produce irritating/toxic fumes/gases.

5.3 Advice for firefighters

Personal protection equipment:

Fire-fighters should wear appropriate protective equipment and self-contained breathing apparatus (SCBA) with a full-face piece operated in positive pressure mode.

Special precautions:

Avoid contact with skin, eyes, hair and clothing. Do not breathe fumes/gas/mists/aerosols/vapors/dusts. Move containers from fire area if safe to do so. Use water spray/fog for cooling fire exposed containers. Avoid unnecessary run-off of extinguishing media which may cause pollution.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures:

Evacuate unnecessary personnel. Ventilate area. Extinguish any sources of ignition. Wear recommended personal protective equipment (see Section 8). Do not get on skin, eyes or on clothing. Avoid breathing mist, vapor, dust, fume and spray. Do not walk through spilled material. Wash thoroughly after handling. Remove contaminated clothing and laundry before reuse.

6.2 Environmental precautions:

Prevent further leakage or spillage if safe to do so. Prevent from reaching drains, sewers and waterways. Discharge into the environment must be avoided.

6.3 Methods and material for containment and cleaning up:

Harmful if swallowed. Do not touch damaged containers or spilled material unless wearing appropriate personal protective clothing. Avoid breathing dust, mist, fumes, vapors or spray. Stop leak if you can do it without risk. Contain and collect spillage and place in suitable container for future disposal. Dispose of in accordance with all applicable regulations (see Section 13).

6.4 Reference to other sections:

For personal protective equipment see Section 8. For disposal see Section 13.

SECTION 7: Handling and storage

7.1 Precautions for safe handling:

Use appropriate personal protective equipment (see Section 8). Use only with adequate ventilation. Avoid breathing mist/vapor/spray/dust. Do not eat, drink, smoke, or use personal products when handling chemical substances. Avoid contact with skin, eyes and clothing. Wash affected areas thoroughly after handling. Keep away from incompatible materials (See Section 10). Keep containers tightly closed when not in use.

Normal precautions for handling chemicals and potentially infectious materials must be observed.

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 5 of 21

Homocysteine R1

7.2 Conditions for safe storage, including any incompatibilities:

Store in cool, dry, well-ventilated location out of direct sunlight. Keep away from food and beverages. Protect from freezing and physical damage. Store away from heat, open flames and other sources of ignition. Keep container tightly sealed. Store away from incompatible materials (See Section 10).
Store between 2-8 °C

7.3 Specific end use(s):

Refer to Section 1 (Recommended Use).

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

Only those substances with limit values have been included below.

Occupational Exposure limit values:

Country (Legal Basis)	Substance	Identifier	Permissible concentration
Denmark	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 3 mg/m ³ (Dust and mist, organic, total dust)
	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	STEL: 6 mg/m ³ (Dust and mist, organic, total dust)
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	STEL: 0.3 mg/m ³
	Ethanol	64-17-5	TWA: 1900 mg/m ³ ([1000 ppm])
	Ethanol	64-17-5	STEL: 3800 mg/m ³ ([2000 ppm])
Finland	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 5 mg/m ³ (Dust and mist, organic, total dust)
	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	15-Minute STEL: 10 mg/m ³ (Dust and mist, organic, total dust)
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Ethanol	64-17-5	15-Minute STEL: 2500 mg/m ³ ([1300 ppm])
	Ethanol	64-17-5	8-Hour TWA-PEL: 1900 mg/m ³ ([1000 ppm])
Germany (TRGS 900)	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 10 mg/m ³ (Dust, inhalable, insoluble particulates)
	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	STEL: 20 mg/m ³ (Dust, inhalable, insoluble particulates)
	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 1.25 mg/m ³ (Dust, respirable, insoluble particulates)
	Sodium azide	26628-22-8	8-Hour TWA: 0.2 mg/m ³
	Sodium azide	26628-22-8	15-Minute STEL: 0.4 mg/m ³

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 6 of 21

Homocysteine R1

Country (Legal Basis)	Substance	Identifier	Permissible concentration
	Ethanol	64-17-5	Limit Value: 380 mg/m ³ ([200 ppm])
Austria	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 10 mg/m ³ (Dust, inhalable)
	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	STEL: 20 mg/m ³ (Dust, inhalable)
	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 5 mg/m ³ (Dust, respirable)
	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	15-Minute STEL: 10 mg/m ³ (Dust, respirable)
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	STEL: 0.3 mg/m ³ (4 x 15 min)
	Ethanol	64-17-5	Ceiling Limit: 3800 mg/m ³ ([2000 ppm])
	Ethanol	64-17-5	MAK TWA: 1900 mg/m ³ ([1000 ppm])
Belgium	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 10 mg/m ³ (Dust, inhalable)
	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 3 mg/m ³ (Dust, respirable)
	Sodium azide	26628-22-8	Ceiling Limit: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Ethanol	64-17-5	8-Hour TWA: 1907 mg/m ³ ([1000 ppm])
France	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 7 mg/m ³ (Dust, inhalable)
	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 3.5 mg/m ³ (Dust, respirable)
	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Ethanol	64-17-5	15-Minute STEL: 9500 mg/m ³ (5000 ppm; [Indicative])
	Ethanol	64-17-5	TWA: 1900 mg/m ³ (1000 ppm; [Indicative])
Hungary	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 10 mg/m ³ (Dust, inhalable)
	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 6 mg/m ³ (Dust, respirable)
	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Ethanol	64-17-5	15-Minute STEL: 3800 mg/m ³
	Ethanol	64-17-5	8-Hour TWA: 1900 mg/m ³
Ireland	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 10 mg/m ³ (Dust, inhalable)

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 7 of 21

Homocysteine R1

Country (Legal Basis)	Substance	Identifier	Permissible concentration
	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 4 mg/m ³ (Dust, respirable)
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Ethanol	64-17-5	15-Minute STEL: 1000 ppm
Poland	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 10 mg/m ³ (Dust, inhalable)
	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Ethanol	64-17-5	8-Hour TWA: 1900 mg/m ³
Spain	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 10 mg/m ³ (Dust, inhalable)
	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 3 mg/m ³ (Dust, respirable)
	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Ethanol	64-17-5	15-Minute STEL: 1910 mg/m ³ ([1000 ppm])
Sweden	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 10 mg/m ³ (Dust, inhalable)
	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 5 mg/m ³ (Dust, respirable)
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	Ceiling Limit: 0.3 mg/m ³
	Ethanol	64-17-5	Level Limit Value: 1000 mg/m ³ ([500 ppm])
	Ethanol	64-17-5	STEL: 1900 mg/m ³ ([1000 ppm])
Germany (MAK)	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 4 mg/m ³ (Dust, inhalable)
	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	8-Hour TWA: 0.3 mg/m ³ (Dust, respirable, granular, bio-resistant)
	2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	15-Minute STEL: 2.4 mg/m ³ (Dust, respirable, granular, bio-resistant)
	Sodium azide	26628-22-8	8-Hour TWA: 0.2 mg/m ³ (inhalable fraction)
	Sodium azide	26628-22-8	15-Minute STEL: 0.4 mg/m ³ (inhalable fraction)
	Ethanol	64-17-5	8-Hour TWA: 380 mg/m ³ ([200 ppm])
Bulgaria	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 8 of 21

Homocysteine R1

Country (Legal Basis)	Substance	Identifier	Permissible concentration
	Sodium azide	26628-22-8	TWA: 0.1 mg/m ³
	Ethanol	64-17-5	TWA: 1000 mg/m ³
Croatia	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Ethanol	64-17-5	8-Hour TWA: 1900 mg/m ³ ([1000 ppm])
Czech Republic	Sodium azide	26628-22-8	Ceiling Limit: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Ethanol	64-17-5	Ceiling Limit: 3000 mg/m ³
	Ethanol	64-17-5	TWA: 1000 mg/m ³
Estonia	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Ethanol	64-17-5	15-Minute STEL: 1900 mg/m ³ ([1000 ppm])
	Ethanol	64-17-5	8-Hour TWA: 1000 mg/m ³ ([500 ppm])
European Union	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³ ([SCOEL])
Greece	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³ (0.1 ppm)
	Sodium azide	26628-22-8	8-Hour TWA: 0.3 mg/m ³ (0.1 ppm)
	Ethanol	64-17-5	8-Hour TWA: 1900 mg/m ³ ([1000 ppm])
Italy	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Ethanol	64-17-5	15-Minute STEL: 1000 ppm
Latvia	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Ethanol	64-17-5	8-Hour TWA: 1000 mg/m ³
Lithuania	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Ethanol	64-17-5	15-Minute STEL: 1900 mg/m ³ ([1000 ppm])
	Ethanol	64-17-5	8-Hour TWA: 1000 mg/m ³ ([500 ppm])
Luxembourg	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 9 of 21

Homocysteine R1

Country (Legal Basis)	Substance	Identifier	Permissible concentration
	Sodium azide	26628-22-8	TWA: 0.1 mg/m ³
Portugal	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Ethanol	64-17-5	8-Hour TWA: 1000 ppm
Romania	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Ethanol	64-17-5	15-Minute STEL: 9500 mg/m ³ ([5000 ppm])
	Ethanol	64-17-5	8-Hour TWA: 1900 mg/m ³ ([1000 ppm])
Slovakia	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Ethanol	64-17-5	15-Minute STEL: 1920 mg/m ³ ([1000 ppm])
	Ethanol	64-17-5	8-Hour TWA: 960 mg/m ³ ([500 ppm])
Slovenia	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Ethanol	64-17-5	8-Hour TWA: 960 mg/m ³ ([500 ppm])
	Ethanol	64-17-5	15-Minute STEL: 1920 mg/m ³ ([1000 ppm])
The Netherlands	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Ethanol	64-17-5	15-Minute STEL: 1900 mg/m ³ ([1000 ppm])
	Ethanol	64-17-5	8-Hour TWA: 260 mg/m ³ ([137 ppm])
United Kingdom	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Ethanol	64-17-5	8-Hour TWA: 1920 mg/m ³ ([1000 ppm])
Cyprus	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
Malta	Sodium azide	26628-22-8	TWA: 0.1 mg/m ³

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 10 of 21

Homocysteine R1

Country (Legal Basis)	Substance	Identifier	Permissible concentration
	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³

Biological limit values:

No biological exposure limits noted for the ingredient(s).

Derived No Effect Level (DNEL):

Ingredient Name: 2-amino-2-(hydroxymethyl)propane-1,3-diol

CAS #: 77-86-1

Workers - Systemic Effects	Acute - Oral	Not determined or not applicable.
	Acute - Inhalation	No hazard identified
	Acute - Dermal	No hazard identified
	Chronic - Oral	Not determined or not applicable.
	Chronic - Inhalation	117.5
	Chronic - Dermal	166.7 mg/kg bw/day
Workers - Local Effects	Acute - Oral	Not determined or not applicable.
	Acute - Inhalation	No hazard identified
	Acute - Dermal	No hazard identified
	Chronic - Oral	Not determined or not applicable.
	Chronic - Inhalation	No hazard identified
	Chronic - Dermal	No hazard identified
General Population - Systemic Effects	Acute - Oral	No hazard identified
	Acute - Inhalation	No hazard identified
	Acute - Dermal	No hazard identified
	Chronic - Oral	8.3 mg/kg bw/day
	Chronic - Inhalation	29 mg/m ³
	Chronic - Dermal	83.3 mg/kg bw/day
General Population - Local Effect	Acute - Oral	Not determined or not applicable.
	Acute - Inhalation	No hazard identified
	Acute - Dermal	No hazard identified
	Chronic - Oral	Not determined or not applicable.
	Chronic - Inhalation	No hazard identified
	Chronic - Dermal	No hazard identified

Ingredient Name: Sodium azide

CAS #: 26628-22-8

Workers - Systemic Effects	Acute - Oral	Not determined or not applicable.
	Acute - Inhalation	Hazard identified but no DNEL available
	Acute - Dermal	Hazard identified but no DNEL available
	Chronic - Oral	Not determined or not applicable.
	Chronic - Inhalation	0.493 mg/m ³
	Chronic - Dermal	0.14 mg/kg bw/day

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 11 of 21

Homocysteine R1

Workers - Local Effects	Acute - Oral	Not determined or not applicable.
	Acute - Inhalation	No hazard identified
	Acute - Dermal	No hazard identified
	Chronic - Oral	Not determined or not applicable.
	Chronic - Inhalation	No hazard identified
	Chronic - Dermal	No hazard identified
General Population - Systemic Effects	Acute - Oral	Hazard identified but no DNEL available
	Acute - Inhalation	Hazard identified but no DNEL available
	Acute - Dermal	Hazard identified but no DNEL available
	Chronic - Oral	0.05 mg/kg bw/day
	Chronic - Inhalation	0.087 mg/m ³
	Chronic - Dermal	0.05 mg/kg bw/day
General Population - Local Effect	Acute - Oral	Not determined or not applicable.
	Acute - Inhalation	No hazard identified
	Acute - Dermal	No hazard identified
	Chronic - Oral	Not determined or not applicable.
	Chronic - Inhalation	No hazard identified
	Chronic - Dermal	No hazard identified

Ingredient Name: Ethanol

CAS #: 64-17-5

Workers - Systemic Effects	Acute - Oral	Not determined or not applicable.
	Acute - Inhalation	No hazard identified
	Acute - Dermal	No hazard identified
	Chronic - Oral	Not determined or not applicable.
	Chronic - Inhalation	380 mg/m ³
	Chronic - Dermal	No hazard identified
Workers - Local Effects	Acute - Oral	Not determined or not applicable.
	Acute - Inhalation	No hazard identified
	Acute - Dermal	No hazard identified
	Chronic - Oral	Not determined or not applicable.
	Chronic - Inhalation	No hazard identified
	Chronic - Dermal	No hazard identified
General Population - Systemic Effects	Acute - Oral	No hazard identified
	Acute - Inhalation	No hazard identified
	Acute - Dermal	No hazard identified
	Chronic - Oral	No hazard identified
	Chronic - Inhalation	114 mg/m ³
	Chronic - Dermal	No hazard identified

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 12 of 21

Homocysteine R1

General Population - Local Effect	Acute - Oral	Not determined or not applicable.
	Acute - Inhalation	No hazard identified
	Acute - Dermal	No hazard identified
	Chronic - Oral	Not determined or not applicable.
	Chronic - Inhalation	No hazard identified
	Chronic - Dermal	No hazard identified

Predicted No Effect Concentration (PNEC):

Ingredient Name: 2-amino-2-(hydroxymethyl)propane-1,3-diol

CAS #: 77-86-1

Environmental Protection Target	PNEC
Fresh water	No hazard identified
Freshwater sediments	No hazard identified
Marine water	No hazard identified
Marine sediments	No hazard identified
Microorganisms in sewage treatment	300 mg/L
Soil (agricultural)	No hazard identified
Air	No hazard identified
Oral (Secondary Poisoning)	No exposure expected

Ingredient Name: Sodium azide

CAS #: 26628-22-8

Environmental Protection Target	PNEC
Fresh water	0.35 µg/L
Freshwater sediments	0.0167 mg/kg sediment dw
Marine water	0.015 µg/L
Marine sediments	0.00072 mg/kg sediment dw
Microorganisms in sewage treatment	30 µg/L
Soil (agricultural)	No hazard identified
Air	No hazard identified
Food chain	No exposure expected

Ingredient Name: Ethanol

CAS #: 64-17-5

Environmental Protection Target	PNEC
Fresh water	0.96 mg/L
Freshwater sediments	3.6 mg/kg sediment dw
Marine water	0.79 mg/L
Marine sediments	2.9 mg/kg sediment dw
Microorganisms in sewage treatment	580 mg/L
Soil (agricultural)	0.63 mg/kg soil dw
Air	No hazard identified
Oral (Secondary Poisoning)	380 mg/kg food

Information on monitoring procedures:

Not determined or not applicable.

8.2 Exposure controls

Appropriate engineering controls:

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 13 of 21

Homocysteine R1

Emergency eye wash stations and safety showers should be available in the immediate vicinity of use or handling. Provide adequate ventilation to maintain the airborne concentrations of vapor, mists, and/or dusts below the applicable workplace exposure limits, while observing recognized national standards (or equivalent).

Personal protection equipment

Eye and face protection:

Safety glasses or goggles. Use eye protection equipment that has been tested and approved by recognized national standards (or equivalent).

Skin and body protection:

Chemical resistant, impervious gloves approved by the appropriate standards. Gloves must be inspected prior to use. Avoid skin contact with used gloves. Appropriate techniques should be used to remove used gloves and contaminated clothing. Personal protective equipment for the body should be selected based on the task being performed and the risks involved and should be approved by a specialist before handling this product. Ensure that all personal protective equipment is approved by recognized national standards (or equivalent).

Respiratory protection:

If engineering controls do not maintain airborne concentrations below the applicable workplace exposure limits, or to an acceptable level (if exposure limits have not been established), a respirator approved by recognized national standards (or equivalent) must be worn.

General hygienic measures:

When handling chemical products, do not eat, drink or smoke. Wash hands after handling, before breaks, and at the end of the workday. Avoid contact with skin, eyes and clothing. Wash contaminated clothing before reuse. Perform routine housekeeping.

Environmental exposure controls:

Emissions from ventilation or work process equipment should be checked to ensure they comply with the requirements of environmental protection legislation.

Product (substance / mixture) related measures to prevent exposure:	Not determined or not applicable.
Instruction measures to prevent exposure:	Not determined or not applicable.
Organisational measures to prevent exposure:	Not determined or not applicable.
Technical measures to prevent exposure:	Not determined or not applicable.

Risk management measures to control exposure:

Not determined or not applicable.

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Physical State	Stable liquids
Color	Clear and colorless to slight yellow in color.
Odor/Odor threshold	Odorless
pH	Not Available
Melting point/freezing point	Not determined or not available.
Initial boiling point/range	Not Available
Flash point (closed cup)	Not determined or not available.
Flammability	Not determined or not available.
Upper flammability/explosive limit	Not determined or not available.
Lower flammability/explosive limit	Not determined or not available.

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 14 of 21

Homocysteine R1

Vapor pressure	Not Available
Relative vapor density	Not determined or not available.
Density	Not determined or not available.
Relative density	Not determined or not available.
Solubilities	Not Available
Partition coefficient (n-octanol/water)	Not determined or not available.
Auto/Self-ignition temperature	Not Available
Decomposition temperature	Not determined or not available.
Kinematic viscosity	Not Available
Particle characteristics	Not determined or not available.

9.2 Other information

9.2.1 Information with regard to physical hazard classes

Explosives	No data available/Not applicable
Flammable gases	No data available/Not applicable
Aerosols	No data available/Not applicable
Oxidizing gases	No data available/Not applicable
Gases under pressure	No data available/Not applicable
Flammable liquids	No data available/Not applicable
Flammable solids	No data available/Not applicable
Self-reactive substances and mixtures	No data available/Not applicable
Pyrophoric liquids	No data available/Not applicable
Pyrophoric solids	No data available/Not applicable
Self-heating substances and mixtures	No data available/Not applicable
Substances and mixtures, which emit flammable gases in contact with water	No data available/Not applicable
Oxidizing liquids	No data available/Not applicable
Oxidizing solids	No data available/Not applicable
Organic peroxides	No data available/Not applicable
Corrosive to metals	No data available/Not applicable
Desensitized explosives	No data available/Not applicable

9.2.2 Other safety characteristics

None.

SECTION 10: Stability and reactivity

10.1 Reactivity:

Not reactive under recommended handling and storage conditions.

10.2 Chemical stability:

Stable under recommended handling and storage conditions.

10.3 Possibility of hazardous reactions:

Hazardous reactions are not anticipated under recommended conditions of handling and storage.

10.4 Conditions to avoid:

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 15 of 21

Homocysteine R1

Extreme heat, open flames, hot surfaces, sparks, ignition sources and incompatible materials.

10.5 Incompatible materials:

Strong oxidizing agents and acids.

10.6 Hazardous decomposition products:

Under normal conditions of storage and use, hazardous decomposition products should not be produced.

SECTION 11: Toxicological information

11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008

Acute toxicity

Assessment:

Harmful if swallowed.

Product data: No data available.

Substance data:

Name	Route	Result
2-amino-2-(hydroxymethyl)propane-1,3-diol	oral	LD50 Rat: >5000 mg/kg
	dermal	LD50 Rat: >5000 mg/kg
Sodium azide	oral	LD50 Rat: 42 mg/kg
	dermal	LD50 Rabbit: 5 mg/kg
	inhalation	LC50 Rat: >0.054 mg/L (4 hr [Dust])
Ethanol	oral	LD50 Rat: 10,470 mg/kg
	inhalation	LC50 Rat: 116.9 mg/L (4 hr [vapor])
	dermal	LD50 Rabbit: 17,100 mg/kg

Skin corrosion/irritation

Assessment: Based on available data, the classification criteria are not met.

Product data:

No data available.

Substance data:

Name	Result
2-amino-2-(hydroxymethyl)propane-1,3-diol	Causes skin irritation.

Serious eye damage/irritation

Assessment: Based on available data, the classification criteria are not met.

Product data:

No data available.

Substance data:

Name	Result
2-amino-2-(hydroxymethyl)propane-1,3-diol	Causes serious eye irritation.
Ethanol	Causes serious eye irritation.

Respiratory or skin sensitization

Assessment: Based on available data, the classification criteria are not met.

Product data:

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 16 of 21

Homocysteine R1

No data available.

Substance data: No data available.

Carcinogenicity

Assessment: Based on available data, the classification criteria are not met.

Product data: No data available.

Substance data: No data available.

International Agency for Research on Cancer (IARC):

Name	Classification
2-amino-2-(hydroxymethyl)propane-1,3-diol	Not Applicable
Sodium azide	Not Applicable
Ethanol	Not Applicable

Germ cell mutagenicity

Assessment: Based on available data, the classification criteria are not met.

Product data: No data available.

Substance data: No data available.

Reproductive Toxicity

Assessment: Based on available data, the classification criteria are not met.

Product data:

No data available.

Substance data: No data available.

Specific target organ toxicity (single exposure)

Assessment: Based on available data, the classification criteria are not met.

Product data:

No data available.

Substance data:

Name	Result
2-amino-2-(hydroxymethyl)propane-1,3-diol	May cause respiratory irritation.

Specific target organ toxicity (repeated exposure)

Assessment: Based on available data, the classification criteria are not met.

Product data:

No data available.

Substance data:

Name	Result
Sodium azide	May cause damage to the brain through prolonged or repeated exposure.

Aspiration toxicity

Assessment: Based on available data, the classification criteria are not met.

Product data:

No data available.

Substance data: No data available.

Information on likely routes of exposure:

No data available.

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 17 of 21

Homocysteine R1

Symptoms related to the physical, chemical and toxicological characteristics:

No data available.

11.2 Information on other hazards

Endocrine disrupting properties:

Substance data: No data available.

Other information:

No data available.

SECTION 12: Ecological information

12.1 Toxicity

Acute (short-term) toxicity

Assessment: Based on available data, the classification criteria are not met.

Product data: No data available.

Substance data:

Name	Result
2-amino-2-(hydroxymethyl)propane-1,3-diol	Aquatic Invertebrates EC50 Daphnia magna: >980 mg/L (48 hr [mobility])
	Aquatic Plants EC50 Pseudokirchneriella subcapitata: 397 mg/L (72 hr [cell number])
	Fish LC50 Freshwater fish: >4000 mg/L (96 hr)
Sodium azide	Fish LC50 Gasterosteus aculeatus: 0.8 mg/L (96 hr)
	Aquatic Plants EC50 Pseudokirchneriella subcapitata: 0.35 mg/L (96 hr [cell number])
Ethanol	Fish LC50 Pimephales promelas: 15,300 mg/L (96 hr)
	Aquatic Invertebrates LC50 Ceriodaphnia dubia: 5012 mg/L (48 hr)
	Aquatic Plants EC50 Chlorella vulgaris: 275 mg/L (72 hr [growth rate])
	Bacteria LC50 Paramecium caudatum: 5,800 mg/L (4 hr)

Chronic (long-term) toxicity

Assessment: Based on available data, the classification criteria are not met.

Product data: No data available.

Substance data:

Name	Result
Ethanol	Aquatic Invertebrates NOEC Daphnia Magna: 9.6 mg/L (10 d [reproduction])
	Fish NOEC Danio rerio: 250 mg/L (5 d)

12.2 Persistence and degradability

Product data: No data available.

Substance data:

Name	Result
2-amino-2-(hydroxymethyl)propane-1,3-diol	Readily biodegradable in water (65.9% biodegradation in 28 days, measured by CO2 evolution).
Sodium azide	Biodegradation studies do not apply to inorganic substances.
Ethanol	The substance is readily biodegradable. 84% degradation measured by O2 consumption, after 20 days.

12.3 Bioaccumulative potential

Product data: No data available.

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 18 of 21

Homocysteine R1

Substance data:

Name	Result
2-amino-2-(hydroxymethyl)propane-1,3-diol	No significant bioaccumulation expected. Log Kow (aquatic species): -2.31
Ethanol	The substance is not expected to bioaccumulate in organisms (estimated BCF: 3).

12.4 Mobility in soil

Product data: No data available.

Substance data:

Name	Result
2-amino-2-(hydroxymethyl)propane-1,3-diol	Mobile in soil with a low potential for adsorption to soil and sediment. Koc at 20 °C: 75
Ethanol	The substance is highly mobile; therefore, adsorption to soil is not expected (log Koc: 0.2).

12.5 Results of PBT and vPvB assessment

Product data:

PBT assessment: This product does not contain any substances that are assessed to be a PBT.

vPvB assessment: This product does not contain any substances that are assessed to be a vPvB.

Substance data:

PBT assessment:

2-amino-2-(hydroxymethyl)propane-1,3-diol	This substance is not PBT.
Sodium azide	PBT assessment does not apply to inorganic substance.
Ethanol	The substance is not PBT.

vPvB assessment:

2-amino-2-(hydroxymethyl)propane-1,3-diol	This substance is not vPvB.
Sodium azide	vPvB assessment does not apply to inorganic substances.
Ethanol	The substance is not vPvB.

12.6 Endocrine disrupting properties

Substance data: No data available.

12.7 Other adverse effects: No data available.

12.8 Hazard to the ozone layer

Assessment: Based on available data, the classification criteria are not met.

Product data: No data available.

Substance data: No data available.

SECTION 13: Disposal considerations

13.1 Waste treatment methods

13.1.1 Product / Packaging disposal:

Dispose of waste product, unused product and contaminated packaging in compliance with federal, state and local regulations. If unsure of the applicable requirements, contact the authorities for information.

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 19 of 21

Homocysteine R1

Waste codes / waste designations according to LoW: Not determined or not available.

13.1.2 Waste treatment-relevant information: Not determined or not available.

13.1.3 Sewage disposal-relevant information: Not determined or not available.

13.1.4 Other disposal recommendations: It is the responsibility of the waste generator to properly characterize all waste materials according to applicable regulatory entities

SECTION 14: Transport information

International Carriage of Dangerous Goods by Road/Rail (ADR/RID)

UN number or ID number	Not regulated
UN proper shipping name	Not regulated
UN transport hazard class(es)	None
Packing group	None
Environmental hazards	None
Special precautions for user	None

International Carriage of Dangerous Goods by Inland Waterways (ADN)

UN number or ID number	Not regulated
UN proper shipping name	Not regulated
UN transport hazard class(es)	None
Packing group	None
Environmental hazards	None
Special precautions for user	None

International Maritime Dangerous Goods (IMDG)

UN number or ID number	Not regulated
UN proper shipping name	Not regulated
UN transport hazard class(es)	None
Packing group	None
Environmental hazards	None
Special precautions for user	None

International Air Transport Association Dangerous Goods Regulations (IATA-DGR)

UN number or ID number	Not regulated
UN proper shipping name	Not regulated
UN transport hazard class(es)	None
Packing group	None
Environmental hazards	None
Special precautions for user	None

Maritime Transport in Bulk according to IMO Instruments

Bulk Name	None
Ship type	None
Pollution category	None

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 20 of 21

Homocysteine R1

IMO hazard class	None
Environmental hazards	None
Material hazardous only in bulk	None
Cargo Group	None

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture.

European regulations

Inventory listing (EINECS): All ingredients are listed or exempt.

REACH SVHC candidate list: None of the ingredients are listed.

REACH SVHC Authorizations: None of the ingredients are listed.

REACH Restriction: None of the ingredients are listed.

Water hazard class (WGK) (Product): Not determined.

Water hazard class (WGK) (Substance):

Ingredient Name	CAS	Class
2-amino-2-(hydroxymethyl)propane-1,3-diol	77-86-1	Water hazard class 1: slightly hazardous to water
Sodium azide	26628-22-8	Water hazard class 2: obviously hazardous to water
Ethanol	64-17-5	Water hazard class 1: slightly hazardous to water

Other regulations

Germany TA Luft: None of the ingredients are listed.

Additional information: Not determined.

15.2 Chemical Safety Assessment

No Chemical Safety Assessment has been carried out for this substance/mixture by the supplier.

SECTION 16: Other information

Abbreviations and Acronyms: None

Classification procedure:

Classification according to Regulation (EC) No. 1272/2008 (CLP)	Method Used
Acute toxicity (oral), category 4	Expert judgement

Summary of classification(s) in section 3:

Skin Irrit. 2	Skin irritation, category 2
STOT SE 3 (RI)	Specific target organ toxicity - single exposure, category 3, respiratory tract irritation
Eye Irrit. 2	Eye Irritation, category 2
Flam. Liq. 2	Flammable liquids, category 2
Acute Tox. 2 (Oral)	Acute toxicity (oral), category 2
Aquatic Acute 1	Acute aquatic hazard, category 1
Acute Tox. 1 (Dermal)	Acute toxicity (dermal), category 1
Acute Tox. 2 (Inh)	Acute toxicity (inhalation), category 2
Aquatic Chronic 1	Chronic aquatic hazard, category 1
STOT RE 2	Specific target organ toxicity - repeated exposure, category 2

Summary of hazard statements in section 3:

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 21 of 21

Homocysteine R1

H315	Causes skin irritation
H335	May cause respiratory irritation
H319	Causes serious eye irritation
H225	Highly flammable liquid and vapour
H300	Fatal if swallowed
H400	Very toxic to aquatic life
H310	Fatal in contact with skin
H330	Fatal if inhaled
H410	Very toxic to aquatic life with long lasting effects
H373	May cause damage to organs (state all organs affected, if known) through prolonged or repeated exposure (state route of exposure if it is conclusively proven that no other routes of exposure cause the hazard)

Summary of EUH Statement(s) in section 3:

EUH032	Contact with acids liberates very toxic gas
--------	---

Disclaimer:

This product has been classified in accordance with EC No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and EC No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878. The information provided in this SDS is correct, to the best of our knowledge, based on information available. The information given is designed only as a guidance for safe handling, use, storage, transportation, and disposal and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials, unless specified in the text. The responsibility to provide a safe workplace remains with the user.

Initial preparation date: 2023-11-20

End of Safety Data Sheet

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 1 of 15

Homocysteine R2

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Product Name: Homocysteine R2

Product code: H7575-R2

1.2 Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses: For in vitro quantitative determination of total homocysteine in human serum and plasma.

Uses advised against: Not determined or not applicable.

Reasons why uses advised against: Not determined or not applicable.

1.3 Details of the manufacturer/supplier of the safety data sheet

Manufacturer:

United States

HORIBA Instruments Incorporated

5449 Research Drive

Canton, MI 48188

734-487-8300

horiba.com

1.4 Emergency telephone number:

United States

HORIBA Instruments Incorporated

1-800-445-9853 (24 hours per day)

France

Organisme de conseil/centre antipoison national

+33 1 45 42 59 59 (24 hours per day)

Portugal

Órgão consultor nacional/Centro Antivenenos

+351 800 250 250 (24 hours per day)

Spain

Centro de información toxicológica/organismo asesor nacional

+34 91 562 04 20 (24 hours per day)

Czech Republic

Národní poradní orgán/toxikologické středisko

+420 224 919 293 (24 hours per day)

Greece

Εθνικό συμβουλευτικό όργανο/Κέντρο Δηλητηριάσεων

+30 210 779 3777 (24 hours per day)

Italy

Organismo ufficiale di consultazione nazionale/Centro antiveleni

+39 06 305 4343 (24 hours per day)

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 2 of 15

Homocysteine R2

Romania

Organism consultativ național/Centru pentru otrăviri
+40 21 3183606 (24 hours per day)

Poland

Krajowa instytucja doradcza/Ośrodek zatruc
+48 22 619 66 54 (24 hours per day)

SECTION 2: Hazard(s) identification

2.1 Classification of the substance or mixture:

Classification according to Regulation (EC) No. 1272/2008 (CLP):

Acute toxicity (oral), category 4

Hazard-determining components of labeling:

Sodium azide

Additional Information: None

2.2 Label elements

Labelling according to Regulation (EC) No 1272/2008 (CLP)

Hazard pictograms:



Signal Word: Warning

Hazard statements:

H302 Harmful if swallowed

Precautionary statements:

P264 Wash skin thoroughly after handling.

P270 Do not eat, drink or smoke when using this product

P301+P312 IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.

P330 Rinse mouth

P501 Dispose of contents into sewer system after diluting with large volumes of water, if in accordance with local regulations.

2.3 Other hazards: None known

SECTION 3: Composition/information on ingredients

3.1 Substance: Not applicable.

3.2 Mixture:

Identification	EU REACH Registration No.	Name	Classification according to Regulation (EC) No. 1272/2008 (CLP)	Weight %

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 3 of 15

Homocysteine R2

CAS number: 26628-22-8 EC number: 247-852-1	-	Sodium azide	Acute Tox. 2 (Oral); H300 Aquatic Acute 1; H400 Acute Tox. 1 (Dermal); H310 Acute Tox. 2 (Inh); H330 Aquatic Chronic 1; H410 STOT RE 2; H373 M-Factor: 1 EUH032	<1
--	---	--------------	--	----

Additional information: None

Full Text of H and EUH statements: See section 16

SECTION 4: First aid measures

4.1 Description of first aid measures

General notes:

Show this Safety Data Sheet to the doctor in attendance. Take precautions to ensure your own safety before attempting rescue. Wear appropriate safety eyewear, gloves, protective clothing and respiratory protection to prevent exposure. See Section 8 of this SDS for personal protective equipment recommendations. Do not use the mouth to mouth method if victim has ingested or inhaled the product. Give artificial respiration with the aid of a pocket mask equipped with a one-way valve or other proper device.

Following inhalation:

If inhaled, remove person to fresh air and place in a position comfortable for breathing. Keep person at rest. If breathing is difficult, administer oxygen. If breathing has stopped, provide artificial respiration. If symptoms develop or persist, seek medical advice/attention.

Following skin contact:

Remove contaminated clothing and shoes. Rinse skin with copious amounts of water [shower] for several minutes. Launder contaminated clothing before reuse. If symptoms develop or persist, seek medical advice/attention.

Following eye contact:

Rinse eyes with plenty of water for several minutes. Remove contact lenses, if present and easy to do so. Protect unexposed eye. If symptoms develop or persist, seek medical advice/attention.

Following ingestion:

If swallowed, DO NOT induce vomiting unless told to do so by a physician or poison control center. Rinse mouth with water. Never give anything by mouth to an unconscious person. If spontaneous vomiting occurs, place on the left side with head down to prevent aspiration of liquid into the lungs. If symptoms develop or persist, seek medical advice/attention.

Self-Protection of the first aider:

Not determined or not available.

4.2 Most important symptoms and effects, both acute and delayed

Acute symptoms and effects:

Acute oral exposure may lead to dizziness, drowsiness, headache, breathing difficulties, nausea, vomiting, abdominal pain, and lowering of consciousness. Adverse effects are dependent on exposure (dose, concentration, contact time).

Delayed symptoms and effects:

Symptoms of exposure may be delayed.

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 4 of 15

Homocysteine R2

4.3 Indication of any immediate medical attention and special treatment needed

Specific treatment:

Not determined or not available.

Notes for the doctor:

Treat symptomatically.

SECTION 5: Firefighting measures

5.1 Extinguishing media

Suitable extinguishing media:

Water mist/fog, carbon dioxide, dry chemical or alcohol resistant foam.

Unsuitable extinguishing media:

Do not use water jet.

5.2 Special hazards arising from the substance or mixture:

Thermal decomposition may produce irritating/toxic fumes/gases.

5.3 Advice for firefighters

Personal protection equipment:

Fire-fighters should wear appropriate protective equipment and self-contained breathing apparatus (SCBA) with a full-face piece operated in positive pressure mode.

Special precautions:

Avoid contact with skin, eyes, hair and clothing. Do not breathe fumes/gas/mists/aerosols/vapors/dusts. Move containers from fire area if safe to do so. Use water spray/fog for cooling fire exposed containers. Avoid unnecessary run-off of extinguishing media which may cause pollution.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures:

Evacuate unnecessary personnel. Ventilate area. Extinguish any sources of ignition. Wear recommended personal protective equipment (see Section 8). Do not get on skin, eyes or on clothing. Avoid breathing mist, vapor, dust, fume and spray. Do not walk through spilled material. Wash thoroughly after handling. Remove contaminated clothing and laundry before reuse.

6.2 Environmental precautions:

Prevent further leakage or spillage if safe to do so. Prevent from reaching drains, sewers and waterways. Discharge into the environment must be avoided.

6.3 Methods and material for containment and cleaning up:

Harmful if swallowed. Do not touch damaged containers or spilled material unless wearing appropriate personal protective clothing. Avoid breathing dust, mist, fumes, vapors or spray. Stop leak if you can do it without risk. Contain and collect spillage and place in suitable container for future disposal. Dispose of in accordance with all applicable regulations (see Section 13).

6.4 Reference to other sections:

For personal protective equipment see Section 8. For disposal see Section 13.

SECTION 7: Handling and storage

7.1 Precautions for safe handling:

Use appropriate personal protective equipment (see Section 8). Use only with adequate ventilation. Avoid breathing mist/vapor/spray/dust. Do not eat, drink, smoke, or use personal products when handling chemical substances. Avoid contact with skin, eyes and clothing. Wash affected areas thoroughly after handling. Keep away from incompatible materials (See Section 10). Keep containers tightly closed when not in use.

Normal precautions for handling chemicals and potentially infectious materials must be observed.

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 5 of 15

Homocysteine R2

7.2 Conditions for safe storage, including any incompatibilities:

Store in cool, dry, well-ventilated location out of direct sunlight. Keep away from food and beverages. Protect from freezing and physical damage. Store away from heat, open flames and other sources of ignition. Keep container tightly sealed. Store away from incompatible materials (See Section 10).
Store between 2-8 °C

7.3 Specific end use(s):

Refer to Section 1 (Recommended Use).

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

Only those substances with limit values have been included below.

Occupational Exposure limit values:

Country (Legal Basis)	Substance	Identifier	Permissible concentration
Austria	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	STEL: 0.3 mg/m ³ (4 x 15 min)
Belgium	Sodium azide	26628-22-8	Ceiling Limit: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
Bulgaria	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	TWA: 0.1 mg/m ³
Croatia	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
Czech Republic	Sodium azide	26628-22-8	Ceiling Limit: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
Estonia	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
European Union	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³ ([SCOEL])
France	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
Germany (MAK)	Sodium azide	26628-22-8	8-Hour TWA: 0.2 mg/m ³ (inhalable fraction)
	Sodium azide	26628-22-8	15-Minute STEL: 0.4 mg/m ³ (inhalable fraction)
Greece	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³ (0.1 ppm)
	Sodium azide	26628-22-8	8-Hour TWA: 0.3 mg/m ³ (0.1 ppm)

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 6 of 15

Homocysteine R2

Country (Legal Basis)	Substance	Identifier	Permissible concentration
Hungary	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
Italy	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
Latvia	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
Lithuania	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
Luxembourg	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	TWA: 0.1 mg/m ³
Poland	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
Portugal	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
Romania	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
Slovakia	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
Slovenia	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
Spain	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
Sweden	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	Ceiling Limit: 0.3 mg/m ³
The Netherlands	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 7 of 15

Homocysteine R2

Country (Legal Basis)	Substance	Identifier	Permissible concentration
United Kingdom	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
Cyprus	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
Malta	Sodium azide	26628-22-8	TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
Denmark	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	STEL: 0.3 mg/m ³
Finland	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³
Germany (TRGS 900)	Sodium azide	26628-22-8	8-Hour TWA: 0.2 mg/m ³
	Sodium azide	26628-22-8	15-Minute STEL: 0.4 mg/m ³
Ireland	Sodium azide	26628-22-8	8-Hour TWA: 0.1 mg/m ³
	Sodium azide	26628-22-8	15-Minute STEL: 0.3 mg/m ³

Biological limit values:

No biological exposure limits noted for the ingredient(s).

Derived No Effect Level (DNEL):

Ingredient Name: Sodium azide

CAS #: 26628-22-8

Workers - Systemic Effects	Acute - Oral	Not determined or not applicable.
	Acute - Inhalation	Hazard identified but no DNEL available
	Acute - Dermal	Hazard identified but no DNEL available
	Chronic - Oral	Not determined or not applicable.
	Chronic - Inhalation	0.493 mg/m ³
	Chronic - Dermal	0.14 mg/kg bw/day
Workers - Local Effects	Acute - Oral	Not determined or not applicable.
	Acute - Inhalation	No hazard identified
	Acute - Dermal	No hazard identified
	Chronic - Oral	Not determined or not applicable.
	Chronic - Inhalation	No hazard identified
	Chronic - Dermal	No hazard identified

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 8 of 15

Homocysteine R2

General Population - Systemic Effects	Acute - Oral	Hazard identified but no DNEL available
	Acute - Inhalation	Hazard identified but no DNEL available
	Acute - Dermal	Hazard identified but no DNEL available
	Chronic - Oral	0.05 mg/kg bw/day
	Chronic - Inhalation	0.087 mg/m ³
	Chronic - Dermal	0.05 mg/kg bw/day
General Population - Local Effect	Acute - Oral	Not determined or not applicable.
	Acute - Inhalation	No hazard identified
	Acute - Dermal	No hazard identified
	Chronic - Oral	Not determined or not applicable.
	Chronic - Inhalation	No hazard identified
Chronic - Dermal	No hazard identified	

Predicted No Effect Concentration (PNEC):

Ingredient Name: Sodium azide

CAS #: 26628-22-8

Environmental Protection Target	PNEC
Fresh water	0.35 µg/L
Freshwater sediments	0.0167 mg/kg sediment dw
Marine water	0.015 µg/L
Marine sediments	0.00072 mg/kg sediment dw
Microorganisms in sewage treatment	30 µg/L
Soil (agricultural)	No hazard identified
Air	No hazard identified
Food chain	No exposure expected

Information on monitoring procedures:

Not determined or not applicable.

8.2 Exposure controls

Appropriate engineering controls:

Emergency eye wash stations and safety showers should be available in the immediate vicinity of use or handling. Provide adequate ventilation to maintain the airborne concentrations of vapor, mists, and/or dusts below the applicable workplace exposure limits, while observing recognized national standards (or equivalent).

Personal protection equipment

Eye and face protection:

Safety glasses or goggles. Use eye protection equipment that has been tested and approved by recognized national standards (or equivalent).

Skin and body protection:

Chemical resistant, impervious gloves approved by the appropriate standards. Gloves must be inspected prior to use. Avoid skin contact with used gloves. Appropriate techniques should be used to remove used gloves and contaminated clothing. Personal protective equipment for the body should be selected based on the task being performed and the risks involved and should be approved by a specialist before handling this product. Ensure that all personal protective equipment is approved by recognized national standards (or equivalent).

Respiratory protection:

If engineering controls do not maintain airborne concentrations below the applicable workplace exposure limits, or to an acceptable level (if exposure limits have not been established), a respirator approved by recognized national standards (or equivalent) must be worn.

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 9 of 15

Homocysteine R2

General hygienic measures:

When handling chemical products, do not eat, drink or smoke. Wash hands after handling, before breaks, and at the end of the workday. Avoid contact with skin, eyes and clothing. Wash contaminated clothing before reuse. Perform routine housekeeping.

Environmental exposure controls:

Emissions from ventilation or work process equipment should be checked to ensure they comply with the requirements of environmental protection legislation.

Product (substance / mixture) related measures to prevent exposure:	Not determined or not applicable.
Instruction measures to prevent exposure:	Not determined or not applicable.
Organisational measures to prevent exposure:	Not determined or not applicable.
Technical measures to prevent exposure:	Not determined or not applicable.

Risk management measures to control exposure:

Not determined or not applicable.

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Physical State	Stable liquids
Color	Clear and colorless to slight yellow in color.
Odor/Odor threshold	Odorless
pH	Not determined or not available.
Melting point/freezing point	Not Available
Initial boiling point/range	Not determined or not available.
Flash point (closed cup)	Not Available
Flammability	Not Available
Upper flammability/explosive limit	Not determined or not available.
Lower flammability/explosive limit	Not determined or not available.
Vapor pressure	Not determined or not available.
Relative vapor density	Not Available
Density	Not determined or not available.
Relative density	Not determined or not available.
Solubilities	Not Available
Partition coefficient (n-octanol/water)	Not determined or not available.
Auto/Self-ignition temperature	Not Available
Decomposition temperature	Not determined or not available.
Kinematic viscosity	Not determined or not available.
Particle characteristics	Not determined or not available.

9.2 Other information

9.2.1 Information with regard to physical hazard classes

Explosives	No data available/Not applicable
Flammable gases	No data available/Not applicable
Aerosols	No data available/Not applicable
Oxidizing gases	No data available/Not applicable

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 10 of 15

Homocysteine R2

Gases under pressure	No data available/Not applicable
Flammable liquids	No data available/Not applicable
Flammable solids	No data available/Not applicable
Self-reactive substances and mixtures	No data available/Not applicable
Pyrophoric liquids	No data available/Not applicable
Pyrophoric solids	No data available/Not applicable
Self-heating substances and mixtures	No data available/Not applicable
Substances and mixtures, which emit flammable gases in contact with water	No data available/Not applicable
Oxidizing liquids	No data available/Not applicable
Oxidizing solids	No data available/Not applicable
Organic peroxides	No data available/Not applicable
Corrosive to metals	No data available/Not applicable
Desensitized explosives	No data available/Not applicable

9.2.2 Other safety characteristics

None.

SECTION 10: Stability and reactivity

10.1 Reactivity:

Not reactive under recommended handling and storage conditions.

10.2 Chemical stability:

Stable under recommended handling and storage conditions.

10.3 Possibility of hazardous reactions:

Hazardous reactions are not anticipated under recommended conditions of handling and storage.

10.4 Conditions to avoid:

Extreme heat, open flames, hot surfaces, sparks, ignition sources and incompatible materials.

10.5 Incompatible materials:

Strong oxidizing agents and acids.

10.6 Hazardous decomposition products:

Under normal conditions of storage and use, hazardous decomposition products should not be produced.

SECTION 11: Toxicological information

11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008

Acute toxicity

Assessment:

Harmful if swallowed.

Product data: No data available.

Substance data:

Name	Route	Result
Sodium azide	oral	LD50 Rat: 42 mg/kg
	dermal	LD50 Rabbit: 5 mg/kg
	inhalation	LC50 Rat: >0.054 mg/L (4 hr [Dust])

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 11 of 15

Homocysteine R2

Skin corrosion/irritation

Assessment: Based on available data, the classification criteria are not met.

Product data:

No data available.

Substance data: No data available.

Serious eye damage/irritation

Assessment: Based on available data, the classification criteria are not met.

Product data:

No data available.

Substance data: No data available.

Respiratory or skin sensitization

Assessment: Based on available data, the classification criteria are not met.

Product data:

No data available.

Substance data: No data available.

Carcinogenicity

Assessment: Based on available data, the classification criteria are not met.

Product data: No data available.

Substance data: No data available.

International Agency for Research on Cancer (IARC):

Name	Classification
Sodium azide	Not Applicable

Germ cell mutagenicity

Assessment: Based on available data, the classification criteria are not met.

Product data: No data available.

Substance data: No data available.

Reproductive Toxicity

Assessment: Based on available data, the classification criteria are not met.

Product data:

No data available.

Substance data: No data available.

Specific target organ toxicity (single exposure)

Assessment: Based on available data, the classification criteria are not met.

Product data:

No data available.

Substance data: No data available.

Specific target organ toxicity (repeated exposure)

Assessment: Based on available data, the classification criteria are not met.

Product data:

No data available.

Substance data:

Name	Result
Sodium azide	May cause damage to the brain through prolonged or repeated exposure.

Aspiration toxicity

Assessment: Based on available data, the classification criteria are not met.

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 12 of 15

Homocysteine R2

Product data:

No data available.

Substance data: No data available.

Information on likely routes of exposure:

No data available.

Symptoms related to the physical, chemical and toxicological characteristics:

No data available.

11.2 Information on other hazards

Endocrine disrupting properties:

Substance data: No data available.

Other information:

No data available.

SECTION 12: Ecological information

12.1 Toxicity

Acute (short-term) toxicity

Assessment: Based on available data, the classification criteria are not met.

Product data: No data available.

Substance data:

Name	Result
Sodium azide	Fish LC50 <i>Gasterosteus aculeatus</i> : 0.8 mg/L (96 hr)
	Aquatic Plants EC50 <i>Pseudokirchneriella subcapitata</i> : 0.35 mg/L (96 hr [cell number])

Chronic (long-term) toxicity

Assessment: Based on available data, the classification criteria are not met.

Product data: No data available.

Substance data: No data available.

12.2 Persistence and degradability

Product data: No data available.

Substance data:

Name	Result
Sodium azide	Biodegradation studies do not apply to inorganic substances.

12.3 Bioaccumulative potential

Product data: No data available.

Substance data: No data available.

12.4 Mobility in soil

Product data: No data available.

Substance data: No data available.

12.5 Results of PBT and vPvB assessment

Product data:

PBT assessment: This product does not contain any substances that are assessed to be a PBT.

vPvB assessment: This product does not contain any substances that are assessed to be a vPvB.

Substance data:

PBT assessment:

Sodium azide	PBT assessment does not apply to inorganic substance.
--------------	---

vPvB assessment:

Sodium azide	vPvB assessment does not apply to inorganic substances.
--------------	---

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 13 of 15

Homocysteine R2

12.6 Endocrine disrupting properties

Substance data: No data available.

12.7 Other adverse effects: No data available.

12.8 Hazard to the ozone layer

Assessment: Based on available data, the classification criteria are not met.

Product data: No data available.

Substance data: No data available.

SECTION 13: Disposal considerations

13.1 Waste treatment methods

13.1.1 Product / Packaging disposal:

Dispose of waste product, unused product and contaminated packaging in compliance with federal, state and local regulations. If unsure of the applicable requirements, contact the authorities for information.

Waste codes / waste designations according to LoW: Not determined or not available.

13.1.2 Waste treatment-relevant information: Not determined or not available.

13.1.3 Sewage disposal-relevant information: Not determined or not available.

13.1.4 Other disposal recommendations: It is the responsibility of the waste generator to properly characterize all waste materials according to applicable regulatory entities

SECTION 14: Transport information

International Carriage of Dangerous Goods by Road/Rail (ADR/RID)

UN number or ID number	Not regulated
UN proper shipping name	Not regulated
UN transport hazard class(es)	None
Packing group	None
Environmental hazards	None
Special precautions for user	None

International Carriage of Dangerous Goods by Inland Waterways (ADN)

UN number or ID number	Not regulated
UN proper shipping name	Not regulated
UN transport hazard class(es)	None
Packing group	None
Environmental hazards	None
Special precautions for user	None

International Maritime Dangerous Goods (IMDG)

UN number or ID number	Not regulated
UN proper shipping name	Not regulated
UN transport hazard class(es)	None
Packing group	None
Environmental hazards	None
Special precautions for user	None

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 14 of 15

Homocysteine R2

International Air Transport Association Dangerous Goods Regulations (IATA-DGR)

UN number or ID number	Not regulated
UN proper shipping name	Not regulated
UN transport hazard class(es)	None
Packing group	None
Environmental hazards	None
Special precautions for user	None

Maritime Transport in Bulk according to IMO Instruments

Bulk Name	None
Ship type	None
Pollution category	None
IMO hazard class	None
Environmental hazards	None
Material hazardous only in bulk	None
Cargo Group	None

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture.

European regulations

Inventory listing (EINECS): All ingredients are listed or exempt.

REACH SVHC candidate list: None of the ingredients are listed.

REACH SVHC Authorizations: None of the ingredients are listed.

REACH Restriction: None of the ingredients are listed.

Water hazard class (WGK) (Product): Not determined.

Water hazard class (WGK) (Substance):

Ingredient Name	CAS	Class
Sodium azide	26628-22-8	Water hazard class 2: obviously hazardous to water

Other regulations

Germany TA Luft: None of the ingredients are listed.

Additional information: Not determined.

15.2 Chemical Safety Assessment

No Chemical Safety Assessment has been carried out for this substance/mixture by the supplier.

SECTION 16: Other information

Abbreviations and Acronyms: None

Classification procedure:

Classification according to Regulation (EC) No. 1272/2008 (CLP)	Method Used
Acute toxicity (oral), category 4	Expert judgement

Summary of classification(s) in section 3:

Acute Tox. 2 (Oral)	Acute toxicity (oral), category 2
Aquatic Acute 1	Acute aquatic hazard, category 1
Acute Tox. 1 (Dermal)	Acute toxicity (dermal), category 1

Safety Data Sheet

According to Regulation (EC) No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and (EC) No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878.

Initial preparation date: 2023-11-20

Page 15 of 15

Homocysteine R2

Acute Tox. 2 (Inh)	Acute toxicity (inhalation), category 2
Aquatic Chronic 1	Chronic aquatic hazard, category 1
STOT RE 2	Specific target organ toxicity - repeated exposure, category 2

Summary of hazard statements in section 3:

H300	Fatal if swallowed
H400	Very toxic to aquatic life
H310	Fatal in contact with skin
H330	Fatal if inhaled
H410	Very toxic to aquatic life with long lasting effects
H373	May cause damage to organs (state all organs affected, if known) through prolonged or repeated exposure (state route of exposure if it is conclusively proven that no other routes of exposure cause the hazard)

Summary of EUH Statement(s) in section 3:

EUH032	Contact with acids liberates very toxic gas
--------	---

Disclaimer:

This product has been classified in accordance with EC No. 1272/2008 (CLP), as amended by Commission Regulation (EU) 2019/521 and Commission Delegated Regulation (EU) 2020/217, and EC No. 1907/2006 (REACH), as amended by Commission Regulation (EU) 2020/878. The information provided in this SDS is correct, to the best of our knowledge, based on information available. The information given is designed only as a guidance for safe handling, use, storage, transportation, and disposal and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials, unless specified in the text. The responsibility to provide a safe workplace remains with the user.

Initial preparation date: 2023-11-20

End of Safety Data Sheet