

# ABX Minocal



| LOT                      |                                                       | CX 490          |        | CAL    |        | <div><div></div><div>(Exp.) 2024-02-05</div><div>(YYYY - MM - DD)</div></div> |        |        |        |                         |         |
|--------------------------|-------------------------------------------------------|-----------------|--------|--------|--------|-------------------------------------------------------------------------------|--------|--------|--------|-------------------------|---------|
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | ABX Lyse        |        |        |        |                                                                               |        |        |        | TOLERANCES<br>TOLERANCE |         |
|                          |                                                       | MICROS          | MICROS | MICROS | MICROS | MICROS                                                                        | PENTRA | PENTRA | PENTRA |                         | PENTRA  |
|                          |                                                       | ES60<br>Care ST |        |        |        | 60 8P<br>60 18P                                                               |        |        |        |                         |         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 9.2             |        |        |        | 9.1                                                                           |        |        |        |                         | ± 0.2   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.49            |        |        |        | 4.41                                                                          |        |        |        |                         | ± 0.06  |
| HB HGB                   | g/dl                                                  | 13.2            |        |        |        | 13.1                                                                          |        |        |        |                         | ± 0.2   |
|                          | g/l                                                   | 132             |        |        |        | 131                                                                           |        |        |        |                         | ± 2     |
|                          | mmol/l                                                | 8.20            |        |        |        | 8.14                                                                          |        |        |        |                         | ± 0.12  |
| HT HCT                   | %                                                     | 37.7            |        |        |        | 37.0                                                                          |        |        |        |                         | ± 1.0   |
|                          | l/l                                                   | 0.377           |        |        |        | 0.370                                                                         |        |        |        |                         | ± 0.010 |
| PLA PLT                  | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 237             |        |        |        | 234                                                                           |        |        |        |                         | ± 10    |
| VMP MPV                  | µm <sup>3</sup> ; fl                                  | 8.6             |        |        |        | 7.7                                                                           |        |        |        |                         | ± 0.5   |

Ref: TEMP-0387 Rev.51 FRONT / RECTO 1300090308

 1300090308  
 Ref: TEMP-0387 Rev.51 FRONT / RECTO