

ABX Minocal



LOT

CX 490

CAL



(Exp.) 2024-02-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60 Care ST				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2				9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.49				4.41					± 0.06
HB HGB	g/dl	13.2				13.1					± 0.2
	g/l	132				131					± 2
	mmol/l	8.20				8.14					± 0.12
HT HCT	%	37.7				37.0					± 1.0
	l/l	0.377				0.370					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	237				234					± 10
VMP MPV	µm ³ ; fl	8.6				7.7					± 0.5

Ref: TEMP-0387 Rev.51 FRONT / RECTO 1300090308