

ABX Minocal



LOT		CX 490		CAL		 (Exp.) 2024-02-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0				9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42				4.41					± 0.06
HB HGB	g/dl	13.1				13.1					± 0.2
	g/l	131				131					± 2
	mmol/l	8.14				8.14					± 0.12
HT HCT	%	37.1				37.7					± 1.0
	l/l	0.371				0.377					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	220				241					± 10
VMP MPV	µm ³ ; fl	8.8				8.9					± 0.5

ABX Minocal



LOT		CX 490		CAL		 (Exp.) 2024-02-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
		MICROS	MICROS								
			Scil Vet ABC Plus ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		8.7								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.41								± 0.06
HB HGB	g/dl		12.8								± 0.2
	g/l		128								± 2
	mmol/l		7.95								± 0.12
HT HCT	%		37.0								± 1.0
	l/l		0.370								± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l		235								± 10
VMP MPV	µm ³ ; fl		8.8								± 0.5

Ref: TEMP-0387 Rev.51 FRONT / RECTO 1300113555