

ABX Minocal

LOT

CX 490

CAL



(Exp.) 2024-02-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0				9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42				4.41					± 0.06
HB HGB	g/dl	13.1				13.1					± 0.2
	g/l	131				131					± 2
	mmol/l	8.14				8.14					± 0.12
HT HCT	%	37.1				37.7					± 1.0
	l/l	0.371				0.377					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	220				241					± 10
VMP MPV	µm ³ ; fl	8.8				8.9					± 0.5

ABX Minocal

LOT

CX 490

Rev 1

CAL



(Exp.) 2024-02-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS	MICROS							
			Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.41							± 0.06
HB HGB	g/dl		12.8							± 0.2
	g/l		128							± 2
	mmol/l		7.95							± 0.12
HT HCT	%		37.0							± 1.0
	l/l		0.370							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l		235							± 10
VMP MPV	µm ³ ; fl		8.8							± 0.5

Ref: TEMP-0387 Rev.51 FRONT / RECTO 1300113555