

ABX Minocal

LOT

CX 490

CAL



(Exp.) 2024-02-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0				9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42				4.41					± 0.06
HB HGB	g/dl	13.1				13.1					± 0.2
	g/l	131				131					± 2
	mmol/l	8.14				8.14					± 0.12
HT HCT	%	37.1				37.7					± 1.0
	l/l	0.371				0.377					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	220				241					± 10
VMP MPV	μm ³ ; fl	8.8				8.9					± 0.5

ABX Minocal

LOT

CX 490

Rev 1

CAL



(Exp.) **2024-02-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS	MICROS							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		Scil Vet ABC Plus ESV 60							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l									± 0.06
HB HGB	g/dl									± 0.2
	g/l									± 2
	mmol/l									± 0.12
HT HCT	%									± 1.0
	l/l									± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l									± 10
VMP MPV	µm ³ ; fl									± 0.5

Ref: TEMP-0387 Rev.51 FRONT / RECTO 1300113555