

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT

MX 445 L
 Rev 1

CONTROL L

(Exp.)

2024-04-05
 (YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST				60 8P - 60 18P ADVIA 60					
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.1				2.1					± 0.4
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	2.18				2.18					± 0.17
HB HGB	g/dl	5.9				5.9					± 0.5
	g/l	59				59					± 5
	mmol/l	3.66				3.66					± 0.31
HT HCT	%	15.5				15.5					± 2.0
	l/l	0.155				0.155					± 0.020
VGM MCV	µm ³ ; fl	71				71					± 4
TGMH MCH	pg	27.1				27.1					± 2.0
	fmol	1.68				1.68					± 0.12
CCMH MCHC	g/dl	38.1				38.1					± 3.0
	g/l	381				381					± 30
	mmol/l	23.67				23.67					± 1.86
IDR RDW	%	16.0				14.5					± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	67				65					± 22
VPM MPV	µm ³ ; fl	8.2				7.5					± 2.0
LYMPHO	#	1.39				1.41					± 0.40
	%	66.0				67.0					± 8.0
MONO	#	0.23				0.20					± 0.20
	%	11.0				9.5					± 6.0
GRANULO	#	0.48				0.49					± 0.40
	%	23.0				23.5					± 7.0

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT

MX 445 N
 Rev 1

CONTROL N

(Exp.)

2024-04-05
 (YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA		PENTRA
		ES60 Care ST				60 8P - 60 18P ADVIA 60					
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.7				7.8					± 0.8
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	4.52				4.56					± 0.26
HB HGB	g/dl	13.4				13.5					± 0.7
	g/l	134				135					± 7
	mmol/l	8.32				8.38					± 0.43
HT HCT	%	37.1				37.4					± 2.7
	l/l	0.371				0.374					± 0.027
VGM MCV	µm ³ ; fl	82				82					± 4
TGMH MCH	pg	29.6				29.6					± 2.0
	fmol	1.84				1.84					± 0.12
CCMH MCHC	g/dl	36.2				36.1					± 3.0
	g/l	362				361					± 30
	mmol/l	22.45				22.42					± 1.86
IDR RDW	%	14.5				13.5					± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	260				260					± 40
VPM MPV	µm ³ ; fl	7.5				7.1					± 2.0
LYMPHO	#	2.81				2.89					± 0.40
	%	36.5				37.0					± 5.0
MONO	#	0.50				0.47					± 0.30
	%	6.5				6.0					± 4.0
GRANULO	#	4.39				4.45					± 0.60
	%	57.0				57.0					± 6.0

Ref: TEMP-0830 Rev.42 FRONT / RECTO 1300104144

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT

MX 445 H
 Rev 1

CONTROL H

(Exp.)

2024-04-05
 (YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA		
		ES60 Care ST				60 8P - 60 18P ADVIA 60						
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	20.6				20.7						± 1.7
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	5.59				5.65						± 0.32
HB HGB	g/dl	17.9				18.1						± 0.9
	g/l	179				181						± 9
	mmol/l	11.12				11.24						± 0.56
HT HCT	%	50.3				50.9						± 3.4
	l/l	0.503				0.509						± 0.034
VGM MCV	µm ³ ; fl	90				90						± 4
TGMH MCH	pg	32.0				32.0						± 2.0
	fmol	1.99				1.99						± 0.12
CCMH MCHC	g/dl	35.6				35.6						± 3.0
	g/l	356				356						± 30
	mmol/l	22.09				22.10						± 1.86
IDR RDW	%	14.0				13.0						± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	512				509						± 65
VPM MPV	µm ³ ; fl	7.7				7.1						± 2.0
LYMPHO	#	3.30				3.11						± 1.00
	%	16.0				15.0						± 5.0
MONO	#	1.65				1.66						± 0.80
	%	8.0				8.0						± 4.0
GRANULO	#	15.66				15.94						± 1.40
	%	76.0				77.0						± 7.0