

ABX Erytrol



LOT

NX 445

CONTROL



(Exp.) 2024-03-05

(YYYY - MM - DD)

Rev 1

| PARAMETRES PARAMETERS | UNITES UNITS | ABX | | | | | | | | | | | |
|--------------------------|---|---------|----------|-------|-------------------------|----------|----------|-------|-------------------------|---------|--------|-------|-------------------------|
| | | CONTROL | | 1 | TOLERANCES TOLERANCE | CONTROL | | 2 | TOLERANCES TOLERANCE | CONTROL | | 3 | TOLERANCES TOLERANCE |
| | | PENTRA | PENTRA | DX120 | | PENTRA | PENTRA | DX120 | | PENTRA | PENTRA | DX120 | |
| | | DX120 | DX NEXUS | | DX NEXUS | DX NEXUS | DX NEXUS | | DX NEXUS | | | | |
| GB WBC | 10 ³ /mm ³ ; 10 ⁹ /l | N/A | 7.4 | | ± 0.8 | N/A | 18.3 | | ± 1.5 | N/A | 2.7 | | ± 0.4 |
| GBC CWBC | 10 ³ /mm ³ ; 10 ⁹ /l | N/A | 7.1 | | ± 0.8 | N/A | 16.7 | | ± 1.4 | N/A | 1.5 | | ± 0.4 |
| ERB | % | N/A | 4.1 | | ± 1.7 | N/A | 9.5 | | ± 3.3 | N/A | 76.9 | | ± 50.0 |
| ERB | # | N/A | 0.291 | | ± 0.100 | N/A | 1.587 | | ± 0.700 | N/A | 1.154 | | ± 0.500 |

1219930076-C

Ref: TEMP-0822 Rev.36

GBC / CWBC

* les valeurs sont données à titre d'indication et doivent être utilisées comme RuO (Utilisation en Recherche Seulement) Indicative expected values, must be used as RuO (Research used only)