

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT	MX 445 L	CONTROL	L	(Exp.)	2024-04-05 (YYYY-MM-DD)
------------	----------	----------------	---	--------	-----------------------------------

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE							
		MICROS															
		EMI CRP															
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l									2.0							± 0.4
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l									2.22							± 0.15
HB HGB	g/dl									6.0							± 0.4
	g/l									60							± 4
	mmol/l									3.73							± 0.25
HT HCT	%									15.5							± 2.0
	l/l									0.155							± 0.020
VGM MCV	µm ³ ; fl									70							± 4
TGMH MCH	pg									27.0							± 2.0
	fmol									1.68							± 0.12
CCMH MCHC	g/dl									38.6							± 3.0
	g/l									386							± 30
	mmol/l									23.98							± 1.86
IDR RDW	%									15.0							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l									57							± 20
VPM MPV	µm ³ ; fl									7.5							± 2.0
LYMPHO	#									1.28							± 0.40
	%									66.5							± 8.0
MONO	#									0.20							± 0.20
	%									7.5							± 6.0
GRANULO	#									0.52							± 0.40
	%									26.0							± 7.0

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MX 445 N

CONTROL N

(Exp.)

2024-04-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		MICROS								
				EMI CRP						
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l			7.9						± 0.8
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l			4.54						± 0.18
HB HGB	g/dl			13.4						± 0.5
	g/l			134						± 5
	mmol/l			8.32						± 0.31
HT HCT	%			36.8						± 2.5
	l/l			0.368						± 0.025
VGM MCV	µm ³ ; fl			81						± 4
TGMH MCH	pg			29.5						± 2.0
	fmol			1.83						± 0.12
CCMH MCHC	g/dl			36.4						± 3.0
	g/l			364						± 30
	mmol/l			22.63						± 1.86
IDR RDW	%			13.5						± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l			255						± 40
VPM MPV	µm ³ ; fl			7.7						± 2.0
LYMPHO	#			3.00						± 0.40
	%			38.0						± 5.0
MONO	#			0.47						± 0.30
	%			6.0						± 4.0
GRANULO	#			4.42						± 0.60
	%			56.0						± 6.0

Ref: TEMP-0830 Rev.42 BACK / VERSO 1300104144

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT		MX 445 H		CONTROL		H		(Exp.)		2024-04-05			
		Rev 1								(YYYY-MM-DD)			
		ABX Lysebio											
PARAMETRES PARAMETERS		UNITES UNITS		MICROS								TOLERANCES TOLERANCE	
				EMI CRP									
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		20.9								± 1.6	
GR	RBC	10 ⁹ /mm ³ ; 10 ¹² /l		5.55								± 0.20	
HB	HGB	g/dl		17.7								± 0.6	
		g/l		177								± 6	
		mmol/l		10.99								± 0.37	
HT	HCT	%		48.8								± 3.0	
		l/l		0.488								± 0.030	
VGM	MCV	µm ³ ; fl		88								± 4	
TGMH	MCH	pg		31.9								± 2.0	
		fmol		1.98								± 0.12	
CCMH	MCHC	g/dl		36.2								± 3.0	
		g/l		362								± 30	
		mmol/l		22.51								± 1.86	
IDR	RDW	%		13.0								± 3.0	
PLAQ.	PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l		490								± 45	
VPM	MPV	µm ³ ; fl		8.1								± 2.0	
LYMPHO		#		3.55								± 1.00	
		%		17.0								± 5.0	
MONO		#		1.67								± 0.80	
		%		8.0								± 4.0	
GRANULO		#		15.68								± 1.40	
		%		75.0								± 7.0	