

# ABX Erytrol



**LOT**

**NX 445**

**CONTROL**



**(Exp.) 2024-03-05**

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL		1	TOLERANCES TOLERANCE	CONTROL		2	TOLERANCES TOLERANCE	CONTROL		3	TOLERANCES TOLERANCE
		PENTRA	PENTRA			PENTRA	PENTRA			PENTRA	PENTRA		
		DX120	DX NEXUS	DX120	DX NEXUS	DX120	DX NEXUS						
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	N/A	7.4		± 0.8	N/A	18.3		± 1.5	N/A	2.7		± 0.4
GBC CWBC	10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	N/A	7.1		± 0.8	N/A	16.7		± 1.4	N/A	1.5		± 0.4
ERB	%	N/A	4.1		± 1.7	N/A	9.5		± 3.3	N/A	76.9		± 50.0
ERB	#	N/A	0.291		± 0.100	N/A	1.587		± 0.700	N/A	1.154		± 0.500

1219930076-C

Ref: TEMP-0822 Rev.36

GBC / CWBC

\* les valeurs sont données à titre d'indication et doivent être utilisées comme RuO (Utilisation en Recherche Seulement) Indicative expected values, must be used as RuO ( Research used only)