

**HORIBA ABX SAS**  
HORIBA ABX SAS  
Parc Euromédecine  
Rue du Caducée - BP 7290  
34184 MONTPELLIER Cedex 4  
FRANCE  
Tél.: 33 (0) 4 67 14 15 16  
Fax.: 33 (0) 4 67 14 15 17

**LOT**

**MC 430**

**CONTROL**



**(Exp.) 2024-04-05**  
(YYYY-MM-DD)

Rev 2

| PARAMETRES<br>PARAMETERS |   | UNITES<br>UNITS |       | ABX Lysebio |         |        |                         |         |         |        |                         |         |         |        |                         |                         |       |
|--------------------------|---|-----------------|-------|-------------|---------|--------|-------------------------|---------|---------|--------|-------------------------|---------|---------|--------|-------------------------|-------------------------|-------|
|                          |   |                 |       | CONTROL     |         |        | 1                       | CONTROL |         |        | 2                       | CONTROL |         |        | 3                       | TOLERANCES<br>TOLERANCE |       |
|                          |   |                 |       | MICROS      |         | PENTRA | TOLERANCES<br>TOLERANCE | MICROS  |         | PENTRA | TOLERANCES<br>TOLERANCE | MICROS  |         | PENTRA | TOLERANCES<br>TOLERANCE |                         |       |
|                          |   |                 |       | EMI CRP     | LC-767G | MS CRP |                         | EMI CRP | LC-767G | MS CRP |                         | EMI CRP | LC-767G | MS CRP |                         |                         |       |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 2.1             | 2.1   | 2.2         | ± 0.4   | 7.9    |                         | 7.6     | 8.2     | ± 1.0  |                         | 20.7    | 20.4    | 21.5   |                         |                         | ± 2.2 |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 2.25            | 2.24  | 2.23        | ± 0.15  | 4.50   | 4.44                    | 4.54    | ± 0.18  | 5.48   | 5.50                    | 5.63    | ± 0.25  |        |                         |                         |       |
| HB HGB                   | g/dl  | 5.9             | 5.9   | 5.7         | ± 0.4   | 13.3   | 13.3                    | 13.5    | ± 0.5   | 17.6   | 17.8                    | 18.1    | ± 0.7   |        |                         |                         |       |
|                          | g/l   | 59              | 59    | 57          | ± 4     | 133    | 133                     | 135     | ± 5     | 176    | 178                     | 181     | ± 7     |        |                         |                         |       |
|                          | mmol/l  | 3.66            | 3.66  | 3.54        | ± 0.25  | 8.26   | 8.26                    | 8.38    | ± 0.31  | 10.93  | 11.05                   | 11.24   | ± 0.44  |        |                         |                         |       |
| HT HCT                   | %   | 16.0            | 15.9  | 15.9        | ± 2.0   | 36.5   | 36.4                    | 36.8    | ± 2.5   | 48.8   | 49.0                    | 49.3    | ± 3.0   |        |                         |                         |       |
|                          | l/l   | 0.160           | 0.159 | 0.159       | ± 0.020 | 0.365  | 0.364                   | 0.368   | ± 0.025 | 0.488  | 0.490                   | 0.493   | ± 0.030 |        |                         |                         |       |
| VGM MCV                  | µm <sup>3</sup> ; fl                                  | 71              | 71    | 71.5        | ± 4     | 81     | 82                      | 81.0    | ± 4     | 89     | 89                      | 87.5    | ± 4     |        |                         |                         |       |
| TGMH MCH                 | pg  | 26.2            | 26.3  | 25.6        | ± 2.4   | 29.6   | 30.0                    | 29.7    | ± 2.4   | 32.1   | 32.4                    | 32.1    | ± 2.8   |        |                         |                         |       |
|                          | fmol  | 1.63            | 1.64  | 1.59        | ± 0.15  | 1.84   | 1.86                    | 1.85    | ± 0.15  | 1.99   | 2.01                    | 2.00    | ± 0.18  |        |                         |                         |       |
| CCMH MCHC                | g/dl  | 36.9            | 37.1  | 35.7        | ± 3.0   | 36.5   | 36.5                    | 36.7    | ± 3.0   | 36.1   | 36.4                    | 36.7    | ± 3.0   |        |                         |                         |       |
|                          | g/l   | 369             | 371   | 357         | ± 30    | 365    | 365                     | 367     | ± 30    | 361    | 364                     | 367     | ± 30    |        |                         |                         |       |
|                          | mmol/l  | 22.94           | 23.04 | 22.20       | ± 1.86  | 22.66  | 22.69                   | 22.80   | ± 1.86  | 22.41  | 22.58                   | 22.82   | ± 1.86  |        |                         |                         |       |
| IDR RDW                  | %   | 15.0            | 14.5  | 16.5        | ± 3.0   | 13.5   | 13.0                    | 13.5    | ± 3.0   | 13.0   | 12.5                    | 12.5    | ± 3.0   |        |                         |                         |       |
| PLAQ. PLTS               | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 59              | 62    | 57          | ± 20    | 250    | 249                     | 256     | ± 40    | 489    | 501                     | 507     | ± 60    |        |                         |                         |       |
| VPM MPV                  | µm <sup>3</sup> ; fl                                  | 7.6             | 8.3   | 8.9         | ± 2.0   | 7.7    | 7.7                     | 8.3     | ± 2.0   | 8.1    | 8.2                     | 8.7     | ± 2.0   |        |                         |                         |       |
| LYMPHO                   | #   | 1.4             | 1.4   | N/A         | ± 0.2   | 2.9    | 2.8                     | N/A     | ± 0.5   | 3.5    | 3.3                     | N/A     | ± 1.1   |        |                         |                         |       |
|                          | %   | 66.5            | 65.5  | N/A         | ± 8.0   | 37.0   | 36.5                    | N/A     | ± 6.0   | 17.0   | 16.0                    | N/A     | ± 5.0   |        |                         |                         |       |
| MONO                     | #   | 0.2             | 0.2   | N/A         | ± 0.2   | 0.5    | 0.5                     | N/A     | ± 0.4   | 2.0    | 2.1                     | N/A     | ± 1.0   |        |                         |                         |       |
|                          | %   | 7.5             | 8.5   | N/A         | ± 7.0   | 6.0    | 6.0                     | N/A     | ± 5.0   | 9.5    | 10.1                    | N/A     | ± 5.0   |        |                         |                         |       |
| GRANULO                  | #   | 0.5             | 0.5   | N/A         | ± 0.2   | 4.5    | 4.4                     | N/A     | ± 0.5   | 15.2   | 15.1                    | N/A     | ± 1.5   |        |                         |                         |       |
|                          | %   | 26.0            | 26.0  | N/A         | ± 7.0   | 57.0   | 57.5                    | N/A     | ± 6.0   | 73.5   | 73.9                    | N/A     | ± 7.0   |        |                         |                         |       |
| CRP                      | mg/dl   | 0.46            | 0.44  | 0.49        | ± 0.40  | 2.43   | 2.39                    | 2.72    | ± 0.80  | 4.42   | 4.33                    | 4.41    | ± 1.50  |        |                         |                         |       |
|                          | mg/l  | 4.6             | 4.4   | 4.9         | ± 4.0   | 24.3   | 23.9                    | 27.2    | ± 8.0   | 44.2   | 43.3                    | 44.1    | ± 15.0  |        |                         |                         |       |

Ref: TEMP-1070 Rev.42 BACK/VERSO 1300092170