

# ABX Minocal



**LOT** CX 472 **CAL** (Exp.) 2022-08-05  
Rev 1 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN								
		H1500 H2500								
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.75</b>								± 0.20
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.54</b>								± 0.06
HB HGB	g/dl	<b>13.5</b>								± 0.2
	g/l	<b>135</b>								± 2
	mmol/l	<b>8.38</b>								± 0.12
HT HCT	%	N/A								N/A
	l/l	N/A								N/A
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>245</b>								± 10
VMP MPV	μm <sup>3</sup> ; fl	N/A								N/A
VGM MCV	μm <sup>3</sup> ; fl	<b>85.0</b>								± 4.0

1300060884  
Ref: TEMP-0387 Rev.48