

ABX Minocal



LOT CX 472 **CAL** (Exp.) 2022-08-05
Rev 1 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN								
		H1500 H2500								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.75								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54								± 0.06
HB HGB	g/dl	13.5								± 0.2
	g/l	135								± 2
	mmol/l	8.38								± 0.12
HT HCT	%	N/A								N/A
	l/l	N/A								N/A
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	245								± 10
VMP MPV	µm ³ ; fl	N/A								N/A
VGM MCV	µm ³ ; fl	85.0								± 4.0

1300060884
Ref: TEMP-0387 Rev.48